**ISOTEST**

USER MANUAL

Title :	<b>ISOTEST</b>		
Document type :	<b>User manual</b>		
Customer reference :		GLC reference:	
Customer :		Number of pages :	74
Version :	C	Date:	07.07.2010

### REVISION HISTORY

Date	Version	Description	Written by
22.12.2009	A	First issue (09/186)	E. Delombaerde
12.02.2010	B	General update (10/037)	S. Chauvin
07.07.2010	C	Issue of software modifications (10/163) – Addition of: <ul style="list-style-type: none"> <li>♦ “Start loading” button</li> <li>♦ “Waiting production phase”</li> <li>♦ a screen to select batch number</li> </ul>	E. Delombaerde

Date	Written by	Reviewed by	Approved by
07.07.2010	E. Delombaerde	S. Chauvin	D. Papin

## CONTENTS

---

	<b>Page</b>
1. IDENTIFICATION.....	5
1.1. ISOLATOR BRAND AND DESCRIPTION.....	5
1.2. VERSION OF THE ISOLATOR AND THE SOFTWARE.....	5
1.3. MANUFACTURER'S NAME AND ADDRESS.....	6
2. COMPLIANCE REPORT .....	7
3. SAFETY NOTES.....	8
3.1. DEFINITIONS AND WARNINGS.....	8
3.2. GENERAL.....	9
4. GENERAL OPERATION AND RANGE OF APPLICATIONS.....	10
5. SAFETY INFORMATION.....	11
6. COMMISSIONING .....	13
7. EQUIPMENT DESCRIPTION .....	14
7.1. OVERVIEW.....	14
7.2. OPERATOR INTERFACE .....	16
7.3. WORK STATION.....	17
7.4. BIODECONTAMINATION AIRLOCK.....	18
7.5. STERITRACE II STERILIZER INTERFACE.....	19
7.6. HANDLING .....	20
8. OPERATION.....	23
8.1. TRANSFER SYSTEM.....	23
8.2. PROCESS .....	25
8.3. STOPPING A CYCLE BY THE USER.....	43
8.4. EMERGENCY AERATION .....	44
8.5. PRINTING.....	45
8.6. PROCESS ALARMS.....	51
8.7. RECIPE.....	52
8.8. CLEANING.....	55

---

9.	INSTRUCTIONS FOR THE TOOLS AND ACCESSORIES .....	56
9.1.	DPTE-BETA BAG™ .....	56
9.2.	TUBING SYSTEM.....	56
10.	OPTION LIST .....	59
10.1.	CONTROL SYSTEM.....	59
10.2.	ACCESSORIES .....	59
10.3.	CONSUMABLES .....	59
11.	DEFINITION OF THE TECHNICAL TERMS .....	60
	APPENDICES .....	61
Appendix 1	"Hydrogen peroxide product sheet" .....	1 page
Appendix 2	"Hydrogen peroxide safety data sheet" .....	11 pages

## 1. IDENTIFICATION


### 1.1. ISOLATOR BRAND AND DESCRIPTION

<b>ISOTEST</b>	Rigid-wall isolator for the performance of sterility tests ( <i>European version</i> )
----------------	---

### 1.2. VERSION OF THE ISOLATOR AND THE SOFTWARE

The isolator is identified by its serial number and by the software version indicated on the operator panel (*see technical manual, configuration chapter*).

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Isolator serial number

Operator panel and PLC program version number

**1.3. MANUFACTURER'S NAME AND ADDRESS**

<b>GETINGE-LA CALHENE</b>	1, rue du Comté de Donegal 41102 Vendôme cedex – France ☎ +33 (0)254 734 747 📠 +33 (0)254 734 710 <a href="http://www.getinge-lacalhene.com">www.getinge-lacalhene.com</a>
---------------------------	--

## 2. COMPLIANCE REPORT

- EMC directive no. 89/336/EEC modified by directives 91/263/EEC, 92/31/EEC, 93/68/EEC, 92/31/EEC.
- Low voltage directive no. 73/23/EEC.
- Test voltage (*standard test*): as per EN 61 010-1, 2001 issue.
- Overvoltage classification II, pollution level 2
- Electrical safety: as per EN 61 010-1, 2001 issue.
- Electromagnetic compatibility:
  - ♦ Emission according to EN 61000-6-4, 2007 issue
  - ♦ Immunity according to EN 61000-6-2, 2005 issue
- Personal protective equipment (*gloves*) directive no. 89/686/EEC.

### 3. SAFETY NOTES

#### 3.1. DEFINITIONS AND WARNINGS

 <b>DANGER</b>	<p>Means that failure to apply the correct safety measures will result in death or serious injury.</p>
 <b>WARNING</b>	<p>Means that failure to apply the correct safety measures may result in death or serious injury.</p>
 <b>CAUTION</b>	<p>Means, when placed next to a warning triangle, that failure to apply the correct safety measures may result in mild injury.</p>
<b>CAUTION</b>	<p>Means, when not placed next to a warning triangle, that failure to apply the correct safety measures may result in material damage.</p>
<b>IMPORTANT</b>	<p>Means that, if the COMMENTS in question are not taken into account, an unwanted result or status may be obtained.</p>
<b>COMMENTS</b>	<p>In these documents, "COMMENTS" are used to draw the reader's attention to important information concerning the product or a specific part of the document.</p>
	<p>Means that the reader should consult the user manual to find out the usage procedures and limitations</p>

### 3.2. GENERAL

 <p><b>WARNING</b></p>	<p>This equipment uses hydrogen peroxide, electricity and compressed air.</p> <p>Failure to observe the recommendations and instructions included in this manual may result in serious injury or considerable material damage. Only qualified personnel completely familiarised with all the safety rules and the installation, operation and maintenance procedures set forth in the various manuals is authorised to work on this equipment. The correct and safe operation of this equipment requires installation, use and maintenance in accordance with good engineering practices.</p>
 <p><b>CAUTION</b></p>	<p>This equipment should only be used for the purposes indicated by its manufacturer. Unauthorised modifications and the use of spare parts and accessories that are not sold or recommended by the equipment manufacturer may result in fires, electric shocks and injuries.</p>
<p><b>IMPORTANT</b></p>	<p>This manual must be kept within reach nearby the equipment, and a copy must be sent to every user. All measurement and testing operations to be performed on the live equipment must comply with the applicable work safety regulations in the country in question. It is advisable to use suitable tools. Before installation and commissioning, please take the time to read the safety instructions and the warnings included in this manual, as well as all the warning labels affixed to the equipment. Please make sure that these warning labels are always readable and replace any missing or damaged labels.</p>
<p><b>Qualified personnel</b></p>	<p>In this manual, "qualified person" refers to a person who is familiar with the installation, assembly, starting-up and use of the equipment, as well as with the risks incurred.</p>
<p><b>Use of the equipment according to its destination</b></p>	<p>The equipment should only be used for the applications specified in the manual, with devices and components recommended and approved by Getinge-La Calhène</p>

#### 4. **GENERAL OPERATION AND RANGE OF APPLICATIONS**

The ISOTEST isolator is a sealed enclosure which allows its user to perform biodecontamination of various surfaces inside the isolator by chemical means in a closed circuit. The unit is designed and equipped so as to be able to be biodecontaminated using hydrogen peroxide, H<sub>2</sub>O<sub>2</sub>.

**APPLICATION:** Rigid-wall isolator for the performance of sterility tests.

The isolator is equipped with a ventilation and filtering module that maintains an overpressure inside the sterile volume while guaranteeing minimum air renewal with HEPA-filtered air. The isolator is furthermore equipped with transfer systems allowing the insertion and removal of sterile equipment and/or products. Handling systems allow the operator to perform operations inside the sterile contained environment from the outside.

These instructions are valid for the ISOTEST (*Rigid-wall stainless-steel isolator with fixed glass window and door*) operating with overpressure and biodecontaminated with hydrogen peroxide, H<sub>2</sub>O<sub>2</sub>. It consists of a stainless-steel biodecontamination airlock with interlocked doors (*inflatable seals*), a working isolator for performing the sterility test (*using a Millipore Equinox pump or equivalent*), and a ventilation/filtering system.

## 5. SAFETY INFORMATION

 <p><b>DANGER</b></p>	<p>Do not use the ISOCYT isolator in an explosive atmosphere. Only use "HYDROCYDE" hydrogen peroxide with a concentration of 35 % as recommended by Getinge-La Calhène</p>
 <p><b>WARNING</b></p>	<p>The ISOCYT isolator makes use of hazardous toxic products. The user must take all necessary steps before the proper authorities to have this new equipment included in their safety analysis study. The technicians who will use these machines must have previously undergone a specific adapted training process.</p> <p>When handling hydrogen peroxide, please consult the product toxicity sheet (<i>see appendix</i>).</p> <p>It is compulsory to wear gloves and goggles when handling hydrogen peroxide.</p> <p>Only persons who have been trained in using the ISOCYT isolator are authorised to perform biodecontamination cycles.</p> <p><b>In the event of a leak of hydrogen peroxide in liquid or gaseous form during a biodecontamination cycle, the operators <b>must press the "EMERGENCY AERATION" button twice in less than 8 seconds and evacuate the hazard area.</b></b></p> <p>An emergency / evacuation procedure must be in place for this purpose.</p> <p>During the biodecontamination phase, the isolators must be connected to an extractor linked directly to the outside of a building to allow aeration.</p> <p>Make sure that the extractors, generally located on the roof, are protected against rain and high winds and sufficiently separated from passage areas and other heating, ventilation and A/C air intakes.</p> <p>The isolator pipes may only be removed after carrying out after the auxiliary isolator aeration phase, taking all necessary precautions in order to compensate for the risk of residual H<sub>2</sub>O<sub>2</sub> in the pipes.</p>

 <p><b>WARNING</b></p>	<p>All maintenance operations must be performed after aerating the circuits that contain hydrogen peroxide in gaseous form and purging the circuits in which the sterilant is in liquid form. It is compulsory to rinse with water when removing a component on the injection circuit.</p>
 <p><b>CAUTION</b></p>	<p>It is advisable to install hydrogen peroxide detectors in the room where the isolator is installed (<i>acceptable exposure threshold: 1 ppm – 1.4 mg/m<sup>3</sup></i>).</p> <p>It is advisable to install a luminous signal to show when a biodecontamination phase is in progress to warn operators who want to access the room (<i>contact available on the terminal block</i>).</p>
 <p><b>CAUTION</b></p>	<p>The outer door of the biodecontamination airlock <b>SHOULD NEVER</b> be used as a support for setting any type of object, even on a temporary basis. This may affect the adjustment of the door, compromising its seal.</p>
 <p><b>DANGER</b></p>	<p>Before opening the inner door, it is advisable to make sure that no objects on the work top will interfere with its opening. It is also imperative for the operator handling the door to make sure that the second operator (<i>using the other station</i>) has been warned and has not placed his hands under the door. When opening, use one hand to hold the handle, making sure to operate the door under optimum conditions. Failure to observe this instruction may result in injuries to the operator due to impact with the door or his hands getting caught between the door and an object placed on the work top.</p>
<p><b>CAUTION</b></p>	<p>When closing one of the two doors (<i>inner door and outer door of the biodecontamination airlock</i>), it is advisable to check that no objects will get caught between the door and the seal. This may cause a leak-tightness fault. The leak test before a biodecontamination cycle may fail in this case, making it impossible to use the device.</p>


## 6. COMMISSIONING

- Set the master switch to position 1.
- Check that the operator panel switches on. The welcome view appears automatically:

## 7. EQUIPMENT DESCRIPTION

### 7.1. OVERVIEW


Number	Designation
1	Work station
2	Sleeves and gloves
3	Quick biodecontamination airlock
4	H <sub>2</sub> O <sub>2</sub> container
5	Red/green LEDs for door opening authorisation
6	Touch-sensitive operator interface

Number	Designation
1	Ventilation / Filtration
2	Control-command
3	Biodecontamination airlock
4	Built-in STERITRACE II hydrogen peroxide sterilizer


## 7.2. OPERATOR INTERFACE

Number	Designation
1	Master power switch
2	<b>EMERGENCY AERATION</b> push-button
3	Touch screen


**7.3. WORK STATION**

Number	Designation
1	Steritest pump
2	DPTE® S transfer system
3	Basket hanging bars
4	Sleeve and glove

 <b>DANGER</b>	Special precautions must be taken when opening the inner door. Consult chapter 5 – Information and Safety.
--	--


**7.4. BIODECONTAMINATION AIRLOCK**

Number	Designation
1	Transfer trolley
2	Biodecontamination airlock door
3	Sealing gaskets
4	Inner door
5	Basket
6	Inspection window

**IMPORTANT** To help when loading the baskets in the transfer trolley, a preparation trolley can be placed in front of the airlock.

### 7.5. STERITRACE II STERILIZER INTERFACE

Number	Designation
1	Needle-positioning lever ( <i>perforation of the cap of the vial</i> )
2	Vial container
3	Hydrogen peroxide vial

## 7.6. HANDLING


Before each sterilization cycle, Getinge-La Calhène recommends checking the integrity of the isolator gloves. The method consists of:

- conducting a test by raising the oxygen level to assess glove performance (*GLT testing device – operator manual NT 3020/1*)

If the result of the test is "not correct", the glove must be replaced.

The procedure for replacing a glove without breach of containment is as follows:

### *Installation of gloves for an isolator under positive pressure*

- Pull the glove (2) to be replaced inside out into the sleeve (1).
- Remove the tightening ring (3).
- Position the new glove (5) parallel to the cuff ring (4) aligning the axis of the thumb with the seal of the sleeve.
- Install the new glove (5) on the cuff ring (4) up from the bead (6) of the glove to be replaced.
- Place the tightening ring (3) around the glove behind the bead (6).
- Position the bead (6) of the glove and the tightening ring (3) on the end of the cuff ring groove (4).
- Remove the glove (2) from the sleeve (1).
- Extract the glove (2).

### Glove storage

- Avoid natural or artificial light.
- Store them in a dry location with temperatures comprised between 16 and 25 ° C.
- Avoid compression stress due to excessive stacking.

### Washing and disinfection

Pay close attention when performing these operations. Wash the gloves with distilled water and then blow dry them at a low temperature.

Disinfection must be performed using a method and products that will not damage the materials (*isopropyl alcohol or another disinfectant product*).


Avoid extended exposure to high temperatures. The use of halogenous products is forbidden, especially in concentrated form (*chloroform, bleach, etc.*)

**Note:** *This list is not exhaustive, and it is always preferable to perform a test on one glove before beginning general application of any new method or product. Glove replacement frequency will be determined according to their use and the frequency with which the isolator is sterilized.*

Manipulation inside the isolator is done through glove-sleeve assemblies. The neoprene gloves are available in different sizes. The sleeves are made from Hypalon® or PVC / PVC-coated polyester.

### CAUTION

Before starting operations, it is compulsory to remove any objects that might damage the glove: watch, rings, etc. Getinge-La Calhène also recommends wearing under-gloves (*for hygiene*) and over-gloves, which must be replaced frequently (*to reduce cross-contamination*).

Number	Designation
1	Gloves
2	R48 O-ring
3	Equipped cuff ring
4	Sleeve
5	R77 O-ring
6	Shoulder ring

***PVC / Divetex sleeve  
Neoprene glove***

***Hypalon® sleeve  
Hypalon® glove***

## 8. OPERATION


### 8.1. TRANSFER SYSTEM

The unit is equipped with two DPTE® ALPHA 190 transfer systems, one of which is fixed to the work top of the work station allowing the connection of a DPTE-DispoBag™ transfer system to ensure safe evacuation of liquid/solid waste. The other one is fixed to the vertical wall for connecting the Tubing system for evacuating tests.


Once the Beta part or parts are connected, the double door or doors can be opened by turning the lever.

Each one of the two work stations (*facing each other*) can be equipped with these DPTE® transfer systems.

It consists of two parts:

Number	Designation
1	DPTE-S® 190 ( <i>Alpha part</i> )
2	Dummy container ( <i>Beta part</i> ) – ( <i>optional</i> )

The following is provided for the connection of a Beta part: DPTE-DispoBag™ or DPTE-Tubing™.


This system is completely interlocked. It allows completely safe connection / disconnection of the Alpha and Beta parts. The door can only be opened using the lever when a Beta part is connected to the Alpha part and, likewise, the Beta part can only be disconnected if the door is closed.

The dummy container is used to facilitate biodecontamination of the seal surface. It can be replaced with a multipurpose container (*supplied as standard*).

### Operating principle of the DPTE® system

#### Description of principle

- 2 Alpha and Beta parts
- 2 opposing lip seals
- 4 overlapping components, interlocked by rotation of the Beta part

## 8.2. PROCESS


After a power failure or when a user logs off, the following view appears:



To access the user view, press the screen.

**User view**

To start a complete cycle, the operator must follow the instructions on the panel with the help of the block diagram.


Number	Designation
1	Operator instructions
2	Block diagram of a cross-section of the machine
3	User logged in
4	Alarm in progress ( <i>appears when an alarm is active</i> ) OR <div style="border: 1px solid black; padding: 2px; display: inline-block; margin-bottom: 5px;"> <span style="background-color: #0056b3; color: white; padding: 2px 5px;">Process values</span> <span style="background-color: #0056b3; color: white; padding: 2px 5px;">Trend</span> <span style="background-color: #0056b3; color: white; padding: 2px 5px;">PID</span> </div> providing access to the process values or the curves or the PID.
5	Buttons to access the various menus. <div style="margin-bottom: 5px;"> <span style="background-color: #0056b3; color: white; padding: 2px 5px; border: 1px solid black;">System</span> Provides access to the system configuration menu (<i>access level &gt;= Maintenance</i>).                     </div> <div style="margin-bottom: 5px;"> <span style="background-color: #0056b3; color: white; padding: 2px 5px; border: 1px solid black;">Recipe</span> Provides access to the recipes (<i>access level &gt;= Validation</i>).                     </div> <div style="margin-bottom: 5px;"> <span style="background-color: #0056b3; color: white; padding: 2px 5px; border: 1px solid black;">Report</span> Provides access to the reports (<i>access level &gt;= Operator</i>).                     </div> <div style="margin-bottom: 5px;"> <span style="background-color: #0056b3; color: white; padding: 2px 5px; border: 1px solid black;">Logout</span> Allows the operator to log out from the system.                     </div> <div style="margin-bottom: 5px;"> <span style="background-color: #0056b3; color: white; padding: 2px 5px; border: 1px solid black;">Alarm</span> Displays the alarm window (<i>access level &gt;= Operator</i>).                     </div>

Number	Designation
6	<p>Isolator control buttons.</p> <p> Halts an isolator cycle in progress (<i>in the case of a biodecontamination, the airlock automatically launches the emergency phase</i>).</p> <p> Launches the isolator cycle.</p> <p> Launches a new batch (<i>in production or waiting production phase</i>).</p> <p> Pressing this button displays the locking conditions.</p>
7	<p>Isolator operation status.</p>
8	<p>Airlock control buttons.</p> <p> Halts an isolator cycle in progress (<i>in the case of a biodecontamination, the airlock automatically launches the emergency phase</i>).</p> <p> Launches the airlock cycle.</p> <p> Pressing this button displays the locking conditions.</p> <p> Allows the user to select a recipe for biodecontamination.</p> <p> Allows to start the loading phase if the automatic transition "Loading → Unloading" is not selected</p>
9	<p>Biodecontamination airlock operation status.</p>
10	<p>Appears if an alarm is active. The  button acknowledges the alarms.</p>
11	<p>Door status LED (<i>red = locked; green = unlocked</i>)</p>
12	<p>Displays the weight and the expiry data of the H<sub>2</sub>O<sub>2</sub> vial.</p>

 <b>CAUTION</b>	<p>The administrator customer is responsible for defining the users and their passwords and for deleting the default users in the following table (see <i>technical manual</i>).</p>
---	--

### 8.2.1. User password

- Operator mode requires the entry of a user password.
- When any of the buttons are pressed, the following screen appears:

- Enter the username (*default "Oper"*) followed by the password (*default "101"*).
- Confirm your entry by pressing the "OK" button.
- The logged-in user can then be seen in the following view:

*Default password table*

	Access level	User	Password
0	Visitor	N/A	N/A
1	Operator	Oper	101
2	Validation	Valid	201
3	Maintenance	Maint	301
9	Administrator	Admin	100

### 8.2.2. Cycle progress


"The goal of steps 1 to 4 is to prepare the isolator during initial commissioning or after a cleaning or maintenance operation with breach of containment (*operation recommended once per month*). The isolator is not sterile at the start of this phase, and is sterile at the end. Steps 5 to 14 describe the sequence of phases during routine production.

#### Step 1:

- Preparing the isolator for biodecontamination.


#### The operations are as follows:

- ♦ Make sure that the inside of the isolator is clean and free from objects (*stainless-steel baskets, vials, equipment, etc.*). Tools designed to remain inside the isolator (*scales, etc.*) can remain there. Objects that are left inside the isolator may risk compromising the quality of the biodecontamination, since the biodecontamination cycle parameters have been validated with the isolator empty.
- ♦ Make sure that the gloves and sleeves are installed correctly.
- ♦ Take a sleeve support and position it in the sleeve and on the RGI 300 by "separating" the spring-mounted ring.
- ♦ "Slip" the cuff ring onto the sleeve support until it comes to a stop against the bead.
- ♦ "Hitch" the sleeve up as far as possible to uncover all its surfaces.
- ♦ Do the same for the other sleeve supports.
- ♦ Perform these steps in the opposite order to remove the sleeve supports.
- ♦ Make sure a dummy container or a multipurpose container is properly connected to each DPTE® system and that the doors of all the DPTE® systems are open.
- ♦ Make sure the outer door of the biodecontamination airlock is closed.


- ♦ Make sure the inner door of the biodecontamination airlock is open.
- ♦ Make sure a vial of HYDROCYDE is inserted in the container of the Steritrac II sterilizer and that the needle-positioning lever is in the bottom position.


**Step 2:** – If the "Start" button is displayed, all the conditions are in place to begin a cycle. Press this button to move on to step 2.

- If the "Not ready" button is displayed, one or more conditions are missing to begin a cycle.

- Press the "Not ready" button and the following window appears.


- All the locking parameters must be set to "YES" in order to begin the cycle.

Yes -> Condition met  
 No -> Condition missing


- Press the enter key to return to the user view.

*Note: Correct the missing condition and then restart the process from step 1.*


**Step 3:** - The cycle starts in test phase

- If the verification of the H<sub>2</sub>O<sub>2</sub> bottle is correct, the following window appears. The following phases then succeed one another automatically (see steps 3 and 4).


- If the cycle takes place without any faults and if the test is declared non-compliant, the following window appears:

- Press "Isolator stop" and return to step 1.

- If the verification of the bottle is not correct, the following window appears.


**Note:** Correct the missing condition and then restart the process from step 1.


- Three possible alarms:
  - ♦ H<sub>2</sub>O<sub>2</sub> quantity for next cycle fault,
  - ♦ H<sub>2</sub>O<sub>2</sub> vial expiry date fault,
  - ♦ H<sub>2</sub>O<sub>2</sub> vial presence fault,
- Replace the H<sub>2</sub>O<sub>2</sub> vial and acknowledge the alarm. A new verification takes place.

**Step 4:** - The test phase is correct.

If the test phase is correct, the biodecontamination cycle starts automatically and the following animation appears:

- Step 5:** – If the biodecontamination phase is completed with no errors, the isolator automatically enters the waiting production phase and the airlock enters the unloading phase.
- Step 6:** – To start a production phase, press  button – the following windows appears:


- Step 7:** – Enter a batch number and press “Start new batch” button – the production phase starts automatically.

**Perform the following operations:**

- ♦ Close all the doors of the DPTE®.
- ♦ Disconnect the multipurpose containers (*or dummy containers*).
- ♦ Connect instead the necessary accessories (*according to the chosen option*): DPTE-Tubing™ or DPTE-DispoBag™.
- ♦ Open the door of the airlock. Note that it will remain open for all the remaining steps of the process.


The following animation appears:

The operator can transfer the load (*in the presence of a load*) into the isolator and then close the door.

**Step 8:** – The isolator door is closed. The airlock is in loading phase.

The airlock is not sterile; the door of the airlock can be opened. The following animation appears:


After opening the door of the airlock, the operator can insert the load.

- Step 9:** – The airlock door is closed. The airlock is in loading phase.

The airlock is not sterile.

The following animation appears:


When the load is inserted, the operator can close the airlock door.

- Step 10:** – The airlock door is closed. The airlock is in loading phase.

The airlock is not sterile. The operator can still open the airlock door (*if opened, the animation in step 6 is displayed*).


The following animation appears (*if there is a transition before starting biodecontamination in manual mode; in the case of "auto" selection, biodecontamination starts automatically after closing the door*):

- The "Start" button is displayed.  
All the conditions required for beginning a cycle are met. Press this button to move on to step 9.
- The operator wants to select another recipe for the biodecontamination phase.  
Press the "Recipe selection" button and the following view appears:


- Select a recipe and then confirm. The following animation appears:


The "Start" button and the selected recipe are displayed. All the conditions required for beginning a cycle are met. Press this button to move on to step 9.

- Step 11:** – The airlock begins to check the H<sub>2</sub>O<sub>2</sub> bottle.


- If the verification of the sensor is not correct, the following window appears.


**Note:** Correct the missing condition and then restart the process from step 7.

- Three possible alarms:
- ♦ H<sub>2</sub>O<sub>2</sub> quantity for next cycle fault,
  - ♦ H<sub>2</sub>O<sub>2</sub> vial expiry date fault,
  - ♦ H<sub>2</sub>O<sub>2</sub> vial presence fault.
- Replace the H<sub>2</sub>O<sub>2</sub> vial and acknowledge the alarm. A new verification takes place.

- If the verification of the bottle is correct, the following window appears. The test phase starts.


- If the test phase is correct, the following phases then succeed one another automatically (see steps 10 and 11).
- If the cycle takes place without any faults and the test is declared non-compliant, the following window appears:


- Check the tightening of the clamps, cable glands and hose attachments.
- Press "stop airlock" and return to step 7.
- If the fault continues, consult the Getinge After-Sales Service Centre for your region.

- Step 12:**
- The test phase is correct.
  - If the leak test phase is correct, the biodecontamination cycle starts and the following animation appears:


- Step 13:**
- If the biodecontamination phase is completed with no errors, the airlock enters the unloading phase. The airlock is now sterile. The following animation appears:


- The operator can transfer the load to the isolator.

- After transferring the load to the isolator, the operator can close the door again.


**Step 14:**

- The isolator door is closed. The airlock is in loading phase. The airlock is not sterile; the door of the airlock can be opened. The following animation appears:


- After opening the door of the airlock, the operator can insert the load.

- Step 15:** – The airlock door is closed. The airlock is in loading phase.  
The airlock is not sterile.  
The following animation appears:


When the load is inserted, the operator can close the airlock door.

- Step 16:** – The airlock door is closed. The airlock is in loading phase.  
The airlock is not sterile. The operator can still open the airlock door (*if opened, the animation in step 12 is displayed*).  
– If the test phase is still valid, the following animation appears:


- If the test phase is not valid, the following animation appears:

- Perform the same steps for subsequent cycles.

### 8.3. STOPPING A CYCLE BY THE USER


- The operator can halt a cycle by pressing the isolator or airlock "stop" button.
- The following view appears:

- If the "YES" button is pressed, the cycle is halted.

 <b>DANGER</b>	Premature halting of the aeration process will allow the outer door to be opened with a residual level of H <sub>2</sub> O <sub>2</sub> in the machine which is higher than the minimum level defined for normal operation. This therefore represents a risk of inhalation for the operator when opening the outer door of the biodecontamination airlock.
--	--

- If the "NO" button is pressed, the cycle continues.
- The following view appears during production:

- If the "YES" button is pressed, the cycle is halted and the isolator is no longer sterile.
- If the "NO" button is pressed, the cycle continues.

#### 8.4. EMERGENCY AERATION

In the event of a problem during the biodecontamination phase, the operator can launch emergency aeration by pressing the red button next to the "EMERGENCY AERATION" touch screen.


## 8.5. PRINTING

If the printing option is selected, a report is automatically generated at the end of each phase.

The data from the preceding cycle are always stored in the memory and can then be printed from the user screen. This is particularly useful for reprinting in the event of a printer error (*no paper, empty ink cartridge, paper jam, etc.*).


### **Procedure for printing a phase report**

- In the user view, press the "report" button and the following view appears:


The screenshot shows the GETINGE user interface. At the top, it displays 'GLC' and the date/time '15/02/2010 11:24:18'. Below this, there are two status indicators: 'AIRLOCK NOT STERILE' and 'ISOLATOR NOT STERILE Stand by'. The main area contains two dropdown menus: 'Airlock report selection' and 'Isolator report selection'. An 'Exit' button is located at the bottom right.

- Select the phase to print (*isolator or airlock*).


The screenshot shows the GETINGE user interface after selecting a phase. It displays 'GLC' and the date/time '15/02/2010 11:24:43'. Below this, there are two status indicators: 'AIRLOCK NOT STERILE' and 'ISOLATOR NOT STERILE Stand by'. The main area displays 'ISOLATOR REPORT' and 'BIO DECONTAMINATION PHASE'. There are three buttons: 'Current cycle', 'Last cycle', and 'Return'. An 'Exit' button is located at the bottom right.

- Select "current cycle" or "last cycle"

- Press the  button to scroll down the menu until the following window appears:

- Press the  button to scroll down the menu until the following window appears:

GLC		GETINGE		
15/02/2010 11:26:49				
<b>AIRLOCK</b>	<b>NOT STERILE</b>	<b>ISOLATOR</b>	<b>NOT STERILE</b>	
		<b>Bio decontamination</b>		
		Min	Max	Ave
<b>STABILIZATION</b>	0	Min		
Airflow (m3/h)		0.0	0.0	0.0
H2O2 injected (g)		0.0		
<b>VAPORIZER AERATION</b>	0	Min		
<b>ISOLATOR AERATION</b>	0	Min		

  
3/3

**Exit**

- Press the  button to print the report.

## GETINGE

Page 1

**BIO DECONTAMINATION**

Isolator Nb	00000
Beginning the	31/12/2002 10:59:59
End	31/12/2002 10:59:59
User	00000000000000000000
Cycle Nb	000000
Batch ID	00000000000000000000

**ISOLATOR**

Pressure set	0000 Pa
	Min    Max    Ave    Stdev
Pressure (Pa)	0000.0    0000.0    0000.0    000.00
Temperature (°C)	000.0    000.0

**CONDITION**

Humidity (%/hr)	00.0
-----------------	------

**PRIMING**

PRIMING	0000 Sec
---------	----------

**INJECTION**

INJECTION	0000 Min
Airflow (m3/h)	000.0    000.0    000.0
H2O2 injected (g)	00.0

**CONDITION 2**

CONDITION 2	0000 Min
-------------	----------

**STABILIZATION**

STABILIZATION	0000 Min
Airflow (m3/h)	000.0    000.0    000.0
H2O2 injected (g)	00.0

**AERATION**

AERATION	0000 Min
----------	----------

**ISOLATOR AERATION**

ISOLATOR AERATION	0000 Min
-------------------	----------


Signature

Perform the same steps for all the subsequent phases.


### Procedure for printing the alarm log

- In the user view, press the "alarm" button and the following view appears:

- Press the "print" button.

- There are two options for printing the alarm log:
  - ♦ All: Prints all the alarms in the memory.
  - ♦ By date: Prints the alarms comprised between the start date and the end date.

– Alarm log

**GETINGE**

Page 1

**Alarms**

**Print the** 31/12/2000 10:59:59 **By** 0000000000000000


Time	Status	Date
12:00:00	KGQ	01/01/1999
	normal <i>kursiv</i>	normal
12:00:00	KGQ	01/01/1999
	normal <b>blinker</b>	normal
12:00:00	KGQ	01/01/1999
	Message Text	
12:00:00	KGQ	01/01/1999
	normal <b>bold</b>	normal
12:00:00	KGQ	01/01/1999
	normal <u>underline</u>	normal

**Alarm status:**

- A alarm appearing.
- (A) D alarm disappearing and acknowledged.

## 8.6. PROCESS ALARMS

- Any error occurring during a cycle is displayed in the alarm view and must be acknowledged.

GLC		GETINGE	
15/02/2010 11:42:29			
AIRLOCK NOT STERILE		ISOLATOR NOT STERILE	
Stand by		Stand by	
Time	Date	Status	
10:51:47	15/02/2010	(A)D	
H2O2 quantity for next cycle fault			
10:49:26	15/02/2010	A	
H2O2 quantity for next cycle fault			
10:18:28	15/02/2010	(A)D	
H2O2 quantity for next cycle fault			
10:11:10	15/02/2010	A	
H2O2 quantity for next cycle fault			
Print		Erase alarms	
Exit			

### Alarm status:

- A alarm appearing (*in red*).
- (A) D alarm disappearing and acknowledged (*in green*).

 Deletes all the alarms in the memory (*access level >= Maintenance*).

**Note:** Please consult the technical manual for the list of alarms and corrective actions.

## 8.7. RECIPE


The parameters are defined and validated by Getinge-La Calhène. Any modification may cause a change in the cycle.

The "recipe menu" can be accessed from the main view with at least "validation" access level.


 <b>WARNING</b>	<p>Do not modify the parameters of the "Isolator biodecontamination", "Airlock biodecontamination", "Isolator emergency" and "Emergency" recipes without first checking the effectiveness of the aeration phase – See technical manual (<i>calibration and verification chapter</i>).</p>
---	---

### Procedure for modifying recipe parameters

- In the user view, press the "recipe" button and the following view appears:


- Select the recipe to be modified.

Test recipe		
Isolator Nb	0	
Inflating pressure	0	Pa
Inflating max time	0	Sec
Stabilization duration	0	Sec
Min test pressure	0	Pa
Pressure drop	0	Pa/min

- Press the value to be modified. The following window opens:

Min: 1      Max: 600

Upper limit value

Lower limit value

Validation after entering value

- Press the  button to print the recipe.

## GETINGE

TEST PHASE RECIPE

Isolator Nb	00000	
Modify the	31/12/2002 10:59:59	
By	000000000000000000	
Inflating pressure	0000	Pa
Inflating max time	000	Sec
Stabilization duration	000	Sec
Min test pressure	0000	Pa
Pressure drop	0000.0	Pa/min
Print the	31/12/2000 10:59:59	
<u>Signature</u>		


Recipe modification date and time

User who modified the recipe

- Perform the same steps for the recipes of the isolator biodecontamination, production, unloading, loading and emergency phases.


**Airlock biodecontamination recipe modification**

- Select the airlock biodecontamination recipe. The following view appears:

Number	Designation
1	Allows the user to select a recipe
2	Number of the selected recipe
3	Name of the parameter to be changed
4	Prints the selected recipe
5	Saves the selected recipe with a new name
6	Saves the modifications to the selected recipe
7	Directly transfers the parameters of the selected recipe to the PLC
8	Deletes the selected recipe

- A confirmation pop-up window appears when saving a recipe.

## 8.8. CLEANING

**It is advisable prior to any cleaning operation on the isolator, to consult the leaflet NT 3015/12.**

<b>CAUTION</b>	Before beginning to use any other cleaning products, it is essential to check their chemical compatibility with the sterilant, the materials used to manufacture the isolator and its various accessories.
----------------	--

## 9. INSTRUCTIONS FOR THE TOOLS AND ACCESSORIES

### 9.1. DPTE-BETA BAG™

Safe waste disposal with DPTE-BetaBag™.

The DPTE® 190 beta is fixed to a mixed bag (20-litre liquid – 100-litre solid). The overall system has been 100 % tested for leaks.

It is available, ready for use and gamma sterilized.

Like all the DPTE® systems, they provide transversal airtight protection between the operator and the environment.


### 9.2. TUBING SYSTEM

The “Tubing” system has been designed to allow dynamic sterile transfer of equipment or products from inside the isolator to the outside, in a semi-continuous way without breach of containment.

The outlet tube uses the DPTE® system and the operating mode is the same as for a standard container.

The tube and DPTE® assembly is sterilized by radiation.

- It consists of two parts:

Number	Designation
1	Tube
2	Welding machine

**9.2.1. Basket**

- The baskets are made of 316 L stainless steel.
- Dimensions (*mm*): 240 x 160 x 370
- Max. admissible load: 2 kg

**9.2.2. Basket ( $\frac{1}{2}$  campaign)**

- The basket is made of 316 L stainless steel.
- Dimensions (*mm*): 240 x 160 x 370
- Max. admissible load: 5 kg


### 9.2.3. Welding machine

The welding machine and its frame are made of 316L stainless steel. The welding machine consists of a double welding line and a cutting system.

This machine makes it possible to isolate the product to be transferred while keeping the containment in the isolator.

The sleeve is cut with precision between the two weld lines.

## 10. OPTION LIST

### 10.1. CONTROL SYSTEM

<input type="checkbox"/>	PLC server traceability
<input type="checkbox"/>	Printer ( <i>with PLC system</i> )
<input type="checkbox"/>	Printer ( <i>without PLC system</i> )

### 10.2. ACCESSORIES

<input type="checkbox"/>	Trolley for mixed waste
<input type="checkbox"/>	4 stainless-steel baskets
<input type="checkbox"/>	H <sub>2</sub> O <sub>2</sub> metering system for the operator room
<input type="checkbox"/>	2 ergonomic seats
<input type="checkbox"/>	1 welding machine ( <i>Tubing™ system only as an option</i> )
<input type="checkbox"/>	2 welding machines ( <i>Tubing™ system only as an option</i> )

### 10.3. CONSUMABLES

<input type="checkbox"/>	DPTE-Beta Bag® 190, mixed waste
<input type="checkbox"/>	Tube ( <i>Tubing system™ only as an option</i> )
<input type="checkbox"/>	10 neoprene gloves, 3/10, Size 7 - Sterile
<input type="checkbox"/>	6 Hypalon sleeves
<input type="checkbox"/>	12 vials of hydrogen peroxide with RFID chips
<input type="checkbox"/>	Cleaning kit

## 11. DEFINITION OF THE TECHNICAL TERMS

<b>Work station (or work isolator)</b>	Sealed volume in which the work is performed.
<b>STERITRACE II</b>	Hydrogen peroxide sterilizer ( $H_2O_2$ ).
<b>Open loop</b>	The sterilizer only takes care of the sterilant intake, the outlet is provided by the isolator extractor.
<b>Closed loop</b>	The sterilizer takes care of the sterilant intake and outlet in the enclosure.
<b>Dust accumulation class</b>	Classification for comparing the airborne contamination levels of one enclosure with another.
<b>Biodecontamination</b>	Sterilization of surfaces using a sterilant in vapour phase.
<b>DPTE® transfer system</b>	Secure transfer system with a double door. The safest method for introducing and eliminating sterile and/or toxic products without breach of containment.
<b>DPTE® ALPHA</b>	Cell clamp / cell door assembly. Fixed part of the DPTE® mounted on the isolator wall.
<b>DPTE® BETA</b>	Container clamp / container door / container body assembly. Mobile part of the DPTE® ( <i>for example a container</i> ).
<b>DPTE-BetaBag™</b>	System for evacuating waste with no risk of contaminating the environment.
<b>Tubing System</b>	The tubing system has been designed to allow dynamic sterile transfer of equipment or products from inside the isolator to the outside, in a semi-continuous way without breach of containment.
<b>Production phase</b>	Work campaign
<b>Aeration phase</b>	Isolator / airlock aeration to eliminate sterilant vapours.
<b>RFID</b>	Radio Frequency Identification. <b>Radio frequency identification</b> is a method for remotely storing and recovering <a href="#">data</a> using markers called RFID tags or RFID transponders. RFID tags are small objects, such as <a href="#">self-adhesive labels</a> , which can be affixed or added to objects or products or even implanted in living organisms. RFID tags consist of an <a href="#">antenna</a> associated with an electronic chip which allows them to receive and reply to radio requests emitted by the transceiver. These electronic chips contain an identifier and, optionally, additional data.

---

# ***APPENDICES***

***(2 appendices)***

---

# GETINGE

## HYDROCYDE – Hydrogen Peroxide 35%

### Specifications

---

	<i>Minimum</i>	<i>Maximum</i>
Hydrogen peroxide, %	35.0	35.7
pH, Apparent	2.4	3.8
Color, APHA		10
Appearance	Clear	

### Properties

---

	<i>Minimum</i>	<i>Maximum</i>
Active Oxygen, %	16.4	23.5
Specific gravity at 20°C, g/mL	9.45	9.98
Boiling point at 760 mm Hg, °C	1.132	1.195
Freezing point, °C	-32	-51
Vapor Pressure at 30°C, mm Hg	18	23
Appearance	Clear liquid with a slightly pungent odor	

### Miscellaneous

---

CAS No.: 7722-84-1

HYDROCYDE meets the specifications for hydrogen peroxide as required by the “Food Chemicals Codex” (FCC), 5th Edition. FCC requirements are as follows:

ASSAY:	Within range stated
LEAD, as Pb:	≤ 4 ppm
RESIDUE ON EVAPORATION:	≤ 60 ppm
PHOSPHATE:	≤ 50 ppm
ACIDITY, as H <sub>2</sub> SO <sub>4</sub> :	< 0.03%
IRON:	< 0.5 ppm
TIN:	< 10 ppm

**BEFORE HANDLING THIS MATERIAL, READ AND UNDERSTAND THE MSDS (MATERIAL SAFETY DATA SHEET) FOR ADDITIONAL INFORMATION ON PERSONAL PROTECTIVE EQUIPMENT AND FOR SAFETY, HEALTH AND ENVIRONMENTAL INFORMATION.**

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#### COMPLETE SOLUTIONS FOR CONTAMINATION PREVENTION

Getinge is the world's leading provider of solutions for effective cleaning, disinfection and sterilization in the healthcare and life science sectors. We are dedicated to helping our customers provide maximum productivity in the most cost-efficient way. We do this by offering well thought through and customized solutions. This means that we are with our customers all the way from architectural planning and education to traceability and support – with complete solutions, long-term commitment and global presence. Getinge – Always with you.

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**GETINGE**  
GETINGE GROUP

[www.getinge.com](http://www.getinge.com)

**GETINGE GROUP** is a leading global provider of equipment and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. **ArjoHuntleigh** focuses on patient handling and hygiene, disinfection, DVT prevention, medical beds, therapeutic surfaces and diagnostics. **GETINGE** provides solutions for infection control and prevention within healthcare and life sciences. **MAQUET** specializes in therapeutic applications, products, solutions and services for OR and ICU.