VESSEL FOR CREAMS – VCR 1002



REACTOR PARA CREMAS – VCR 1002

Fabrication N° 721021-70-1
N° fabricación:

MELT TANK – MET 0400

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Fabrication N°: 721021-70-3
N° fabricación:

Instructions and Maintenance Manual 
Manual de Instrucciones y Mantenimiento 

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Model VCR 1002 / MELT TANK – Model MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

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TABLE OF CONTENTS

1 CONTROL OF CHANGES.....	5
2. SAFETY AND USER INFORMATION	6
2.1 IMPORTANT INFORMATION	6
2.2 GRAPHIC SYMBOLS	6
2.3 SECURITY WARNINGS.....	7
2.4 SAFETY DEVICES.....	8
2.4.1 Interior Overpressure Protection	8
2.4.2 Excess of pressure protection in dimple jacket	8
2.4.3 Emergency Stop	8
2.4.4 Warnings and Alarms	9
3 SUPPLIED EQUIPMENT DESCRIPTION.....	10
4 TECHNICAL CHARACTERISTICS OF CREAM’S VESSEL (Mod. VCR)	10
4.1 GENERAL DESCRIPTION OF VESSEL	10
4.2 COMPONENTS AND ACCESSORIES	14
4.2.1 Manhole and light.....	15
4.2.2 Anchor mixer.....	16
4.2.3 Homogenizer	16
4.2.4 Heating and cooling system	19
4.2.5 Recirculation pipe.....	19
4.2.6 Weighing and dosage system.....	20
4.2.7 Cleaning system.....	20
4.2.8 Venting system	21
4.2.9 Safety elements	22
5 TECHNICAL CHARACTERISTICS OF MELT TANK (Mod. MET)	23
5.1 GENERAL DESCRIPTION OF TANK	23
5.2 COMPONENTS AND ACCESSORIES	26
5.2.1 Mixer system	27
5.2.2 Emptying system	27
5.2.3 Heating system	28
5.2.4 Safety elements	28
6 TECHNICAL CHARACTERISTICS OF SKID FOR TEMPERATURE CONTROL.....	29
6.1 GENERAL DESCRIPTION OF SKID	29
6.2 SKID COMPONENTS AND ACCESSORIES.....	31
7 ELECTRICAL AND CONTROL SYSTEM.....	34
7.1 ELECTRIC SWITCHBOARD PANEL	34
7.2 ELEMENTS INSIDE PANEL	36
7.3 CONTROL PANEL.....	37
7.4 COMMUNICATION SYSTEM.....	38
8 INSTALLATION.....	39
8.1 SECURITY WARNINGS.....	39
8.2 TRANSPORT	39
8.3 EQUIPMENT INSTALLATION.....	39
8.3.1 Previous verification	39
8.3.2 Equipment installation	40
8.3.3 Installation of Weigh Modules in VCR 1002 vessel	40

Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
MDI VCR (v1.0)	

- 9 START-UP AND OPERATION 42**
 - 9.1 START OF OPERATION..... 42**
 - 9.2 RECIRCULATION OPERATIONS 42**
 - 9.3 PRELIMINARY SETTINGS 47**
 - 9.3.1 Setting the drain times.....47
 - 9.3.2 Spill setting47
 - 9.3.3 Setting the maximum weight oscillation48
 - 9.4 EXTERNAL MEMORY 48**
 - 9.5 LIGHT OPERATION..... 48**
 - 9.6 TOUCH-SCREEN OPERATION 49**
 - 9.6.1 Program structure49
 - 9.6.2 OVERVIEW screen (Process / Recipe)50
 - 9.6.3 RECIPE screen53
 - 9.6.4 ALARMS screen55
 - 9.6.5 CHARTS screen.....56
 - 9.6.6 PROGRAM SETTINGS screen.....57
 - 9.6.7 HMI SETTINGS screen.....62
 - 9.6.8 USER SETTINGS screen63
 - 9.7 COMPLIANCE WITH 21 CFR, PART 11 64**
 - 9.7.1 Introduction.....64
 - 9.7.2 Identification and authentication of users64
 - 9.7.3 Automatic log-off.....64
 - 9.7.4 Actions register65
 - 9.7.5 Audit Trail evaluation66
- 10 MAINTENANCE AND CLEANING 68**
 - 10.1 BASIC MAINTENANCE 68**
 - 10.2 MOTORS 69**
 - 10.3 GEAR UNITS..... 70**
 - 10.4 TIGHTNESS SEALS..... 70**
 - 10.5 DRY MECHANICAL SEAL..... 71**
 - 10.6 EQUIPMENT CLEANING 72**
 - 10.7 SPARE PARTS LIST – VCR 1002 74**
 - 10.8 SPARE PARTS LIST – MET 0400 76**

Fabrication Nº:	721021-70-1 / 721021-70-3	Client:	UNIPHARMA
Equipment:	VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400		
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL			MDI VCR (v1.0)

1 CONTROL OF CHANGES

Detail of the document versions to date:

Revision	Date	Description / Reason of Change	Drafted by
00	09/2014	Initial Version	AIRPLAN, S.A.
01	11/2015	Update of section 5 (P&ID)	AIRPLAN, S.A.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

2. SAFETY AND USER INFORMATION

2.1 IMPORTANT INFORMATION

These assembly and maintenance instructions apply exclusively to the equipment detailed in section 3.1, manufactured by AIRPLAN as per the specifications of Unipharma.

This instruction manual describes the equipment’s characteristics and provides instructions for its proper operation, in addition to containing safety notes to be taken into consideration whenever the equipment is operating.

The personnel responsible for assembly, start-up and maintenance must be familiar with these instructions before beginning.

The specifications in this manual must be followed. AIRPLAN SA takes no responsibility for any damage that may be caused by not following these instructions, or by undue use of the equipment.

2.2 GRAPHIC SYMBOLS

This manual contains a series of icons, whose purpose is to highlight the importance of certain explanations and they are described below:

ICON	DESCRIPTIÓN
	<u>Warning: electrical danger</u> Warns that the described operation implies a risk of physical harm if not carried out respecting safety regulations.
	<u>Warning: General danger</u> The described operation implies a risk of physical harm if not carried out respecting safety regulations.
	<u>Information:</u> Highly important information.
	<u>Warning:</u> Provides information that must be heeded to avoid harming the equipment or personnel.

	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

2.3 SECURITY WARNINGS

To maintain the safety of the operator and to avoid possible damage to the equipment or the installation, this manual must be thoroughly read before any operation.

These instructions for use and maintenance include important safety warnings for proper use of the equipment and the installation in general. Following these warnings helps to avoid danger, to reduce maintenance and repair costs and equipment inactivity.

The user must take the appropriate accident prevention measures, respecting that which is stipulated by regulations in force in the country where the equipment is installed. The dangers related to the nature of the products that may be handled should also be taken into account.

Assembly, start-up, use and maintenance of the equipment must be exclusively entrusted to qualified and trained personnel.

This instruction manual must be kept whole and perfectly legible and all operators in charge of using the equipment or responsible for its maintenance must know where the manual is located so as to be able to consult it at all times.

Before using the equipment, ensure that the network's voltage matches the voltage shown on its plaque.

Electrical installation must heed the Low Voltage Regulation in force, and may only be handled by personnel expressly qualified to this end.

	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

2.4 SAFETY DEVICES

2.4.1 Interior Overpressure Protection

The vessel has a rupture disc to protect against an overpressure inside, which is in contact with internal atmosphere. In case of overpressure due to an abnormal operation of the system that exceeds the nominal pressure of rupture disc, it breaks down and depressurizes the vessel.

Interior pressure in vessel is always visible by a pressure gauge located in the top of the vessel.

To ensure the properly operation of the rupture disc, the pressure inside the vessel must not exceed the maximum pressure PS (3 bar). So, it may be necessary the installation of pressure reducer elements which will ensure that the supply pressure is lower than the design pressure.

The user must ensure that these elements are always in perfect operating conditions.

2.4.2 Excess of pressure protection in dimple Jacket

The vessels' double jacket is equipped with internal overpressure protective equipment consisting of a safety valve and a pressure gauge.

In case of overpressure exceeding the safety valve's nominal pressure, the valve would open, thereby relieving the excess pressure inside the vessel's double jacket.

To ensure the properly operation of the security valve the pressure inside the vessel must not exceed the maximum pressure PS (3 bar). So, it may be necessary the installation of pressure reducer elements which will ensure that the supply pressure is lower than the design pressure.

The user must ensure that these elements are always in perfect operating conditions.

IMPORTANT

The relief outlet must be leaded to a safe area out of the reach of any person.

2.4.3 Emergency Stop

The system has an EMERGENCY STOP BUTTON located in the control panel, which stops immediately the equipment when it is pressed.

After an emergency stop, it is necessary to reset the button by turning it on its axis.

Fabrication Nº:	721021-70-1 / 721021-70-3	Client:	UNIPHARMA
Equipment:	VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400		
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL			MDI VCR (v1.0)

2.4.4 Warnings and Alarms

System can generate alarms and warnings:

- **ALARMS:** Risk situation by a breakdown or an unintended equipment operation. In this case, equipment stops immediately and the user is informed about the situation. For return equipment to its normal operation it is required an action performed by user.
- **WARNING:** Non critical situation generated because a parameter has exceeded the limits configured by user (temperature, level,...) The user is informed about situation but the equipment is not stopped.

ALARMS	DESCRIPTION
Emergency stop	It is generated when the emergency stop button in control panel is pressed.
VCR Mixing switch Fuser Mixing switch	It is generated when the VCR or Fuser mixing switch is disconnected.
Homogenizer switch	It is generated when the homogenizer switch is disconnected.
Vacuum pump switch	It is generated when the vacuum pump switch is disconnected.
Temperature VCR / Fuser	It is generated when a fault in temperature sensor operation is detected or it is disconnected.
Inverter failure	It is generated when a fault in homogenizer/VCR mixer/Fuser mixer inverter is detected.
Modbus failure	It is generated when a fault in Modbus communication is detected.
Low Pressure Air	It is generated when inlet pressure in skid is below 0'5 bar.
Weight controller failure	It is generated when a fault in weight controller operation is detected or it is disconnected..
Valve VGPTXX does not open	It is generated when a fault in valve VGPTXX is detected.
Records failure	It is generated when it is no possible to save the generated records (no SD or it is full).

WARNING	DESCRIPTION
The vacuum pump will not start because it is emptying	It is generated when user presses 'VACUUM' button while VCR vessel is emptying.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

3 SUPPLIED EQUIPMENT DESCRIPTION

The following equipment/elements are supplied:

Equipment	Model	Units	Fabrication Nº
Vessel for creams (1000 l)	VCR 1002	1	721021-70-1
Temperature control Skid	TSK SW 09	1	721021-70-2
Melt tank (400 l)	MET 0400	1	721021-70-3
Electrical switchboard	-	1	721021-07
Control panel	-	1	-

4 TECHNICAL CHARACTERISTICS OF CREAM'S VESSEL (MOD. VCR)

4.1 GENERAL DESCRIPTION OF VESSEL

The supplied equipment is designed to creams fabrication, according to the norms ASME BPE Standards (2012) and GMP Guide (2002), and it's equipped with a anchor mixer and an homogenizer which allows to perform the following operations:

- Product homogenization
- Work as recirculation pump during manufacturing process
- Work as cleaning pump during CIP

The vessel consists of a vertical cylindrical body, with welded conical bottom head and torispherical upper head and with a dimple jacket for product cooling and heating in temperature control process.

For protection and safety of people, the equipment has a triple jacket of rock wool and external stainless steel sheet, polished and welded to the vessel's body.

The equipment is suitable to work under vacuum conditions (-1 Kg./cm²) and until 3 kg/cm² inside the vessel and jacket (has U certification according to ASME-BPVC).

In the upper head, there are the following nozzles and accessories:

- Manhole with sight glass and light
- Nozzles for spray balls of CIP/SIP system (3)
- Nozzles for products inlet (purified water, additives)
- Nozzle for rupture disc
- Nozzle for vacuum
- Nozzle for pressure gauge
- Nozzle for venting filter with diaphragm valve
- Nozzle for upper recirculation inlet
- Spare nozzles (3)

In the vessel's body, there are the following connections:

- Double jacket Inlet connection
- Nozzle for lower recirculation inlet
- Nozzle for double jacket accessories (pressure gauge and safety valve)

In the bottom head, there are the following nozzles and accessories:

- Double jacket outlet connection
- Homogenizer connection
- PT-100 Temperature probe

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

The vessel has 3 legs for its support, made from stainless steel. Each leg is equipped with a loading cell for the product dosage control inside the vessel.

Materials:

The main material used in vessel manufacturing has been stainless Steel, 316L quality in all the parts in contact with product and 304L quality in the rest.

The vessel surfaces which are in contact with product are mirror polished with a roughness $Ra < 0,5\mu m$. The external surface are grain polished 180 s/DIN 69100 with a roughness $Ra < 1,2\mu m$.

The seals which are in contact with product are made of EPDM, VITON or SILICONE, according to the FDA standards.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
		MDI VCR (v1.0)

Datasheet / Ficha de datos técnicos
Vessel for creams – VCR Model
Reactor de cremas – Modelo VCR

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¹ The same electric switchboard panel controls VCR 1002 and MET 0400.

4.2 COMPONENTS AND ACCESSORIES

VCR 1002 Model


Figure 01. Sections and top view of VCR vessel with elements and nozzles identification

NOZZLES AND CONNECTIONS

A	Mixer connection flange
B1, B2, B3	Spray Ball
C	Pressure Gauge
D	Rupture disc
E	Venting Filter
F	Spare nozzle
G	Purified Water inlet
H	Vacuum
I	Additives inlet
J1	Upper double jacket inlet
K2	Lower double jacket outlet
L	Temperature detector PT-100
M	Spare Nozzle
N	Lower Recirculation inlet
O	Safety Valve
P	Pressure Gauge
Q	Spare nozzle
R	Upper Recirculation inlet
S	Homogenizer connection flange

ELEMENTS

01	Mixer
02	Spray ball accessory
03	Manhole with sight glass
04	Pressure Gauge
05	Rupture disc
06	Venting filter + valve
07	Temperature detector PT-100
08	Safety valve
09	Pressure Gauge
10	Load Cells

NOZZLES AND CONNECTIONS SUMMARY

Description	Units	Type	Dimensions
Spray ball	3	Clamp ASME BPE	2 ½"
Pressure Gauge	1	Clamp ASME BPE	1 ½"
Rupture Disc	1	Clamp ASME BPE	1 ½"
Venting Filter	1	Clamp ASME BPE	1 ½"
Spare Nozzle	1	Clamp ASME BPE	1 ½"
Purified water inlet	1	Clamp ASME BPE	1 ½"
Vacuum	1	Clamp ASME BPE	1 ½"
Additives inlet	1	Clamp ASME BPE	1 ½"
Upper double jacket inlet	1	W.N. (R.F.) ANSI 816.5	1 ¼"
Lower double jacket outlet	1	W.N. (R.F.) ANSI 816.5	1 ¼"
Temperature Probe PT-100	1	Clamp ASME BPE	1"
Spare Nozzle	1	Clamp ASME BPE	2"
Lower recirculation inlet	1	Clamp ASME BPE	2"
Safety valve	1	BSP	1"
Pressure Gauge	1	BSP	½"
Spare Nozzle	1	Clamp ASME BPE	1 ½"
Upper recirculation inlet	1	Clamp ASME BPE	2"

4.2.1 Manhole and light
MANHOLE WITH SIGHTGLASS AND LIGHT

	Manhole Brand: FRATELLI LAVEGGI Model: D13 Diameter: DN 400 Material: Stainless steel 316L
	Light Brand: Max Müller AG Model: HL125H Light bulb: 24 V, 50 W Units: 1


Equipment can be inspected from the sight glass or manhole. If for any reason the user decides to access to the interior of equipment through the manhole, it is **VERY IMPORTANT** to consider the following:

- Disconnect the equipment and take all necessary measures in order to prevent an accidental start of equipment while people are inside.
- Customer should know the products that use and, in case of toxic or hazardous products, should proceed to take appropriate safety measures.

4.2.2 Anchor mixer

Vessel has an anchor mixer with sanitary Teflon wipers, which facilitates the product homogenization and heat transmission.

BLADE MIXER		
	Manufacturer: BACHILLER	
	Type: Anchor	
	Motor	
	Brand: BALDOR	
	Model: SUPER E	
	Phases: 3	
	Frequency [Hz]: 60	
	Tension [V]: 460	
	Consumption [A]: 4,1	
	Power [kW]: 2,2 kW (3 HP)	
	rpm: 1.755	
	Direction of rotation: Clockwise	
Gear unit		
Brand: NORD		
Model: SK 42125		
Reduction ratio: 18,80		
Sealing:		
Type: Mechanical seal		
Brand: BURGMANN CARTEX SNO		
Materials: Graphite (FDA)/SiC (FDA)/EPDM(FDA)/C4/1.4571		

4.2.3 Homogenizer

The lower stirring system consists of a homogenizer (micronizer), which allows performing the following operations:

- Product homogenization
- Work as recirculation pump during manufacturing process
- Work as cleaning pump during CIP

HOMOGENIZER		
	Manufacturer: IKA WERKE	
	Model: DBI 2000/10	
	Motor:	
	Brand: WEG	
	Model: W22 IE2	
	Nº of Phases: 3	
	Frequency [Hz]: 60	
	Power [kW]: 22	
	Connection [V]: 460	
	Consumption [A]: 34.3	
	rpm: 3.550	

Recirculation without shearing should only last a short period of time since during said period the liquid will not circulate during the underlying dispersion stage and the temperature can rise a great deal in a short amount of time.

	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

Mechanical seal

The body of the pump has a dual-decompression action mechanical seal in the form of a cartridge.

The following materials have been used:

Stationary ring	Rotating ring	O-ring	Springs	Box components	Stationary ring atmosphere	Rotating ring atmosphere
Silicon carbide (Q1)	Silicon carbide (Q1)	FPM – FDA (V16)	Hastelloy C-4 (M)	CrNiMo-Steel 1.4571 (G)	Food-grade graphite (B)	Silicon carbide (Q1)
Alternative:	Food-grade graphite (B)	EPDM (E) or perfluorinated rubber (K)	-	-	-	-

The homogeniser cannot work without cooling the mechanical locks.

Mechanical seal cooling system

The mechanical seal is lubricated and cooled by means of a barrier pressure installation that operates according to the thermosyphon principle and is made up of two pressure hoses and an air tank with the following fittings:

- Pressure gauge
- Gas filling valve
- Inspection hole
- Purge cock
- Barrier fluid filling valve
- Cooling connections
- Level depth switch
- Contact pressure gauge

BARRIER PRESSURE INSTALLATION	
	<p>Manufacturer: AERNI AG</p> <p>Model: AR 139</p> <p>Maximum service pressure: 16 bar</p> <p>Tank volume: 3L</p> <p>Useful volume: 2.5L</p> <p>Service temperature: 200°C</p> <p>Ventilation surface: 0.1m²</p> <p>Cooling coil content: 0.25L</p>

The barrier fluid is water containing 5% of propylene glycol.

The pressure in the seal circuit should exceed by approximately 2bar the working pressure in the pump's chamber.

For information on the barrier fluid's characteristics and the filling conditions, please refer to the homogeniser's manual in the "Manufacturer Manuals" section.

Safety Devices of Barrier System

Barrier system has the following safety devices to ensure the proper operation of the mechanical seal's cooling system.

There is a temperature probe that allows the control of barrier system's temperature and, in case of detecting an excessive temperature, enables automatically the passage of cold water through the barrier system or, in case of detecting maximum temperature, the equipment is stopped automatically.

TEMPERATURE PROBE

Brand: WIKA
 Type: Pt100; Class B
 Model: TR10-C+T19
 Measuring range: 0...+140 °C
 Probe diameter: Ø 3 mm
 Material: Stainless Steel AISI 316L
 Connection: G ½"

To ensure that the barrier system has enough pressure for its operation, there is a pressure gauge with switch contacts, for fluid pressure visualization in barrier system. In case of detecting a pressure under the set pressure, the equipment is stopped automatically.

PRESSURE GAUGE WITH SWITCH CONTACTS

Brand: WIKA
 Model: PGS23.100
 Measuring range: 0 - 10 bar
 Material: Stainless Steel AISI 316L
 Connection: Rosca G ½"

4.2.4 Heating and cooling system

Heating and cooling equipment process is performed by the pass of steam or chilled water through the vessel's dimple jacket

At the bottom of the vessel there is a thermocouple 4-20 mA that sends the temperature parameters to PLC and allows the user to visualize and control the temperature inside the vessel. Temperature control is achieved by opening or closing the steam valve located in the skid.

RESISTANCE TEMPERATURE DETECTOR

Brand:	WIKA
Type:	Pt100; Class B
Model:	TR10-B+T19
Measuring range:	0...+140 °C
Probe diameter:	Ø 6 mm
Material:	Stainless Steel AISI 316L
Connection:	Clamp 1 ½"

All the valves involved in the control and regulation of the cooling and heating process has been integrated in a stainless steel structure (Skid) that is described extensively in the chapter 5 of this manual.

4.2.5 Recirculation pipe

The vessel has a 2" AISI 361L stainless steel pipe with an internal roughness Ra < 0.6µm for product recirculation.

The pipe is equipped with 4 manually operated arc valves that permit the following operations according to their position:

- Recirculating product with shearing which is loaded through the top of the reactor.
- Recirculating product without shearing which is loaded through the top of the reactor.
- Recirculating product with shearing which is loaded through the bottom of the reactor.
- Recirculating product without shearing which is loaded through the bottom of the reactor.
- Recirculation with shearing for CIP.
- Recirculation without shearing for CIP.
- Emptying the reactor.

VALVES

Make:	LIAGLIAG
Model:	ARC
Actuation:	Manual
Quantity:	4
Connection:	2" ASME-BPE clamp

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

It is essential for manually operated valves to be in the right position according to the desired product circulation operation. Section 9.2 of this Manual details the positions of the valves according to each operation.

The tube also has a nozzle with a clamp ferrule for product inlet under vacuum (if desired).

4.2.6 Weighing and dosage system

Each leg of vessel is equipped with a weigh module that sends the information of the vessel's weigh to PLC, from where it is possible to control the dosage of auxiliary products.

WEIGH MODULES	
	Brand: METTLER TOLEDO Model: SBW505 SS Material: Stainless Steel 304 Load Cells: 0745A Cell capacity: 2.200 Kg Units: 3

4.2.7 Cleaning system

There are three cleaning balls on the vessel's cover whose design allows liquid cleaning products (water, detergent) to be introduced in the reactor such that the entire reactor inner wall is covered so as to wash away any product remains stuck thereto.

CLEANING BALLS	
	Manufacturer: LECHLER Model: 5MI.074.1Y.TF07 Material: AISI 316L Quantity: 3

4.2.8 Venting system

The equipment has a sterile and permanently hydrophobic venting filter with a stainless steel housing, located at the top of the vessel, with a porosity of 0,2 µm, that ensures its sterility and prevents microbiological contamination inside the vessel.

VENTIG FILTER WITH HOUSING

Brand:	MILLIPORE
Cartridge model:	Aervent (CTGR75S01)
Pore dimension:	0,20 µm
Housing model:	S3000 (H1UV7S5G1HPP)
Housing material:	316L Stainless steel
Connections:	Clamp ASME-BPE 1 ½"
Units:	1

In the venting filter nozzle there is a diaphragm valve which allows the venting operation when opened and the pressurize operation and vacuum when closed.

VALVE

Brand:	Bürkert
Type:	Diaphragm manual
Model:	3233
Materials:	Body: 316L stainless steel Diaphragm: EPDM
Connections:	Clamp ASME BPE 1 ½"
Dimensions:	DN 40

The filtering groups must not be used with liquids incompatible with the construction material, including cleaning agents. Incompatible materials are those that chemically attack, soften, create cracks from pressure or adversely affect construction materials in any way. Users must verify the compatibility of the liquids used before using them.

The user is responsible for verifying the actual operation conditions to ensure that the group is compatible with the application and local safety regulations.

4.2.9 Safety elements

The equipment has the following safety elements:

PRESSURE GAUGE (VESSEL)

Located in the upper head of vessel allows the visualization of the vacuum and the pressure inside the vessel.

Brand: WIKA
 Model: 233.50
 Diameter: 100 mm
 Material: Inox AISI 316L
 Connection to process: Diaphragm seal
 Range: -1 a +5 bar
 Connection: ASME BPE Clamp 1 ½"

RUPTURE DISC

Located in the upper head of vessel in order to protect against vessel overpressure.

Brand: FIKE
 Model: SR-H
 Rupture pressure: 3 bar (g)
 Rupture temperature: 136 °C
 Material: AISI 316L + EPDM seal
 Connection: Clamp 1 ½"

PRESSURE GAUGE (SKID)

Located in the connection pipe between skid and vessel it shows the jacket pressure.

Brand: WIKA
 Model: 233.50.100
 Range: 0 a 4 bar
 Connection: ½" GAS nut
 Material: AISI 316L

SAFETY VALVE

Equipment has a safety valve. In case of overpressure exceeding the safety valve's nominal pressure, the valve will open, thereby relieving the excess pressure.

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	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

5 TECHNICAL CHARACTERISTICS OF MELT TANK (MOD. MET)

5.1 GENERAL DESCRIPTION OF TANK

The equipment is designed for melting and warming the fat phase that will be subsequently introduced in cream's vessel by vacuum suction.

The melt tank consists of a vertical cylindrical body that has a welded torispherical lower head and flat cover and is provided with a double jacket for product heating through the passage of steam. For personnel protection and safety, it has an external insulating sheet made of rock wool and a polished stainless-steel sheet, welded to the body.

The tank has been designed to work at atmospheric pressure. The double jacket has been designed to work at 3 kg/cm² pressure (has U certification according to ASME-BPVC).

The following components are located on the upper head of the tank:

- Anchor mixer
- Steam regulation valve
- Pressure gauge (in contact with double jacket)
- Safety valve (in contact with double jacket)

The following components are located on the lower head of the tank:

- Butterfly valve for emptying
- PT100 temperature sensor

The tank sits above the ground on 4 stainless-steel equidistant tubular legs with wheels that allows its mobility.

The different parts of the tank are made in AISI 316L stainless steel for components that are in contact with product and AISI 304 stainless steel for the rest.

The tank surfaces that will come in contact with the product have been polished to a mirror-like finish with a roughness Ra < 0.8µm, whereas the other surfaces have been polished to a 180-grain satin finish as per DIN 69100 with a roughness Ra < 1.2µm.

The gaskets that will come in contact with the product are made of EPDM, VITON or silicone suitable for the pharmaceutical and cosmetic industries.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
		MDI VCR (v1.0)

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Tanque fusor - Modelo MET 0400

PRINCIPALES CARACTERÍSTICAS <i>MAIN CHARACTERISTICS</i>	MET 0400
Fondo superior <i>Upper head</i>	Elíptico <i>Elliptical</i>
Fondo inferior <i>Bottom head</i>	Elíptico <i>Elliptical</i>
Cuerpo <i>Body</i>	Cilíndrico <i>Cylindrical</i>
Altura total <i>Total Height</i>	≈ 1638 mm 64,57''
Diámetro interior <i>Interior Diameter</i>	≈ Ø 750 mm 29,53''
Diámetro exterior <i>Exterior Diameter</i>	≈ Ø 868 mm 34,17''
Peso en vacío <i>Unladen Weight</i>	≈ 530 kg
Capacidad Total en reactor <i>Total Capacity in vessel</i>	400 l gal
Capacidad Útil en reactor <i>Useful Capacity in vessel</i>	400 l gal
Capacidad Total en doble cámara <i>Total Capacity in jacket</i>	100 l gal
Temperatura Máxima (interior y doble cámara) <i>Maximum Temperature (interior and jacket)</i>	145°C 293°F
Presión máxima interior reactor <i>Maximum interior Pressure</i>	-
Presión máxima Doble Cámara <i>Maximum jacket Pressure</i>	3bar(g)
MATERIALES Y ACABADOS <i>MATERIALS AND FINISHES</i>	
Cuerpo / Body	
Material / <i>Material</i>	AISI 316 L
Espesor / <i>Thickness</i>	5 mm / 0,2''
Fondo superior / Upper Head	
Material / <i>Material</i>	AISI 316 L
Espesor / <i>Thickness</i>	6 mm / 0,24''
Fondo inferior / Bottom Head	
Material / <i>Material</i>	AISI 316 L
Espesor / <i>Thickness</i>	6 mm / 0,24''
Doble cámara / Dimple jacket	
Material / <i>Material</i>	AISI 316 L
Espesor / <i>Thickness</i>	3 mm / 0,11''
Triple cámara / Triple jacket	
Material / <i>Material</i>	AISI 304 L
Espesor / <i>Thickness</i>	3 mm / 0,11''
Aislamiento / <i>Insulation</i>	Lana de roca / Rock Wool
Acabados / Finishes	
En contacto con el product / <i>In contact with product</i>	Mirror polished, Ra < 0,5 µm
Exterior / <i>Exterior</i>	Grain Polished 180 s/DIN69100, Ra < 1,2

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

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² The same electric switchboard panel controls VCR 1002 and MET 0400.

5.2 COMPONENTS AND ACCESSORIES

MET Model


Figure 1. Side view and top view of the tank

CONNECTIONS AND NOZZLES

A	Mixer connection flange
B	Empty nozzle
C	Upper double jacket inlet
D	Lower double jacket outlet
E	PT-100 Temperature sensor
F	Spare nozzle
G	Safety valve connection
H	Pressure gauge connection

COMPONENTS

01	Empty valve
02	Mixer
03	Temperature sensor PT-100
04	Safety valve
05	Pressure Gauge

NOZZLES AND CONNECTIONS SUMMARY

Description	Units	Type	Dimensions
Empty nozzle	1	Clamp ASME BPE	2 ½"
Upper double jacket inlet	1	W.N. (R.F.) ANSI 816.5	1"
Lower double jacket outlet	1	W.N. (R.F.) ANSI 816.5	1"
PT100 Temperature sensor	1	Clamp ASME BPE	1"
Spare Nozzle	2	Clamp ASME BPE	2"
Safety valve connection	1	BSP	1"
Pressure Gauge Connection	1	BSP	½"

5.2.1 Mixer system

The tank has an anchor mixer for melted product homogenization, equipped with a frequency inverter for adjusting the mixer's rotation speed.

ANCHOR MIXER

	Manufacturer:	BACHILLER
	Type:	Anchor
	Motor	
	Brand:	NORD
	Model:	SK 90 SH4
	Nº of phases:	3
	Frequency [Hz]:	60
	Connection [V]:	230 / 460
	Consumption [A]:	4,30/ 2,13
	Power [kW]:	1,1
rpm:	1.740	
Rotation sense:	Clockwise/Anticlockwise	
Gear		
Brand:	NORD	
Model:	SK 9012	
Reduction Ratio:	24,53	
Sealing		
Type	Double seal with chamber	
Materials	PTFE / EPDM	

5.2.2 Emptying system

The tank is provided with a butterfly valve, manually operated, that allows emptying the tank.

BUTTERFLY VALVE

	Brand:	INOXPA
	Model:	V48A1
	Diameter:	2"
	Materials:	In contact with product: AISI 316L
		Gaskets: EPDM

5.2.3 Heating system

The product is heated by the passage of steam through the double jacket. The lower portion of the vessel has a 4-20mA temperature sensor that allows regulating the temperature acting on the pneumatic steam inlet valve.

TEMPERATURE SENSOR

Manufacturer: WIKA
 Type: PT100; Class B
 Sensor model: TR10-B
 Transmitter model: T19
 Measurement range: 0...+140°C
 Rod material: AISI 316L stainless steel
 Connection: Clamp 1½"

5.2.4 Safety elements

The tank has the following safety elements:

PRESSURE GAUGE

The tank is provided with a pressure gauge at the top of the double chamber to continuously measure the pressure therein.

Manufacturer: WIKA
 Model: 233.50.100
 Range: 0 to 4bar
 Connection: ½" GAS thread
 Material: AISI 316L

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The safety valve is in contact with the atmosphere inside the double chamber. In case of overpressure exceeding the safety valve's nominal pressure, the valve will open, thereby relieving the excess pressure.

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6 TECHNICAL CHARACTERISTICS OF SKID FOR TEMPERATURE CONTROL

6.1 GENERAL DESCRIPTION OF SKID

With the VCR vessel and MET tank is also supplied a skid which integrates all the valves, pneumatic actuators and accessories needed for the cooling and heating process.

Skid consists in a stainless steel structure of 1050x520x1682 mm (Length x Width x Height) and has the following connections:

- SKID Chilled water outlet
- SKID Chilled water inlet
- SKID Steam inlet
- SKID Drain outlet
- MET Steam outlet
- VCR Steam outlet/SKID Chilled Water Inlet
- SKID Condensates outlet
- MET Steam inlet
- VCR Steam Inlet/SKID Chilled Water Outlet

Valves control is performed by a PLC and a touch screen and for the temperature regulation in heating process there are pneumatic regulation valves (4-20 mA).

The skid has also a compressed air inlet in order to ensure the properly drainage of the piping.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
		MDI VCR (v1.0)

Datasheet

Temperature control Skid - TSK SW Model

PRINCIPALES CARACTERÍSTICAS MAIN CHARACTERISTICS	TSK SW09
Longitud Total <i>Total Length</i>	1286 mm <i>50,63"</i>
Altura Total <i>Total Height</i>	1682 mm <i>66,22"</i>
Ancho Total <i>Total Width</i>	520 mm <i>20,47"</i>
MATERIALES Y ACABADOS MATERIALS AND FINISHES	
Estructura / <i>Structure</i>	AISI 304 L
Conducciones / <i>Pipes</i>	AISI 304 L
Aislamiento conducciones / <i>Pipe insulation</i>	Black Elastomeric Insulation
CHILLED WATER	
Presión Entrada Agua Skid <i>Water Inlet Pressure at skid</i>	2 bar(g)
Caudal (agua de chiller 8°C) <i>Flow (chilled water 46'4 °F)</i>	3 m ³ /h
STEAM	
Presión vapor Entrada skid <i>Steam Inlet Pressure at skid</i>	3 bar(g)
Caudal ³ <i>Flow</i>	120 kg/h (VCR) 50 kg/h (MET)
COMPRESSED AIR	
Presión Máxima <i>Maximum Pressure</i>	0,5 bar

³ Recommended steam flow for heating using 2 sections of the dimple jacket.

6.2 SKID COMPONENTS AND ACCESSORIES


Figure 02. Skid interior view with elements identification

Id	Tipo de elemento	Descripción	DN
VGPP01-70	Angle-seat control valve Inox + Control Head	Steam regulation	32
VGPP02-70	Angle-seat control valve Inox + Control Head	Steam regulation	25
VGPT01-70	Angle-seat valve Inox NC + Pneumatic Actuator	Chilled water inlet	32
VGPT05-70	Angle-seat valve Inox NC + Pneumatic Actuator	Chilled water outlet	32
VGPT13-70	Angle-seat valve Inox NC + Pneumatic Actuator	Condensates outlet	25
VGPT14-70	Angle-seat valve Inox NC + Pneumatic Actuator	Condensates outlet	15
VGPT17-70	Angle-seat valve Inox NC + Pneumatic Actuator	Drain outlet	32
VGPT21-70	Angle-seat valve Inox NC + Pneumatic Actuator	Compressed Air inlet	15
CT01	Thermodynamic steam trap Inox with filter PMA 63 bar	-	25
CT02	Thermodynamic steam trap Inox with filter PMA 63 bar	-	15
AT01	Automatic Air eliminator (inox)	-	15

Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
MDI VCR (v1.0)	

ANGLE-SEAT VALVE PNEUMATICALLY ACTUATED

	<p>It corresponds to valves with reference VGPT01-70, VGPT05-70, VGPT13-70, VGPT14-70, VGPT17-70 y VGPT21-70:</p> <p>Manufacturer: Bürkert Model: 2000 DN: 25 for VGPT13-70 15 for VGPT21-70 32 for VGTP01-70, VGPT05-70 y VGPT17-70</p>
---	---

	<p>All valves have an electrical position feedback</p> <p>Manufacturer: Bürkert Model: 1062</p>
---	---

STEAM REGULATION VALVE

	<p>It corresponds with the valves with references VGPP01-70 and VGPP02-70:</p> <p>Manufacturer: Bürkert Model: 8802-YC (valve 2702 + control head 8694) Body material: AISI 316 L Diameter: DN 32 for VGPP01-70 DN 25 for VGPP02-70 Pressure: From 5,5 to 7 bar Electrical supply: 24 VDC; 4-20 mA</p>
--	--

AIR ELIMINATOR

	<p>It corresponds to the element with reference AT01:</p> <p>Manufacturer: SPIRAX-SARCO Model: AE36A Material: AISI 316I DN: 15</p>
---	---

THERMODYNAMIC STEAM TRAP


	<p>It corresponds to the elemnt with references CT01 and CT02:</p> <p>Manufacturer: SPIRAX-SARCO Model: TD42H Material: Stainless steel DN: 25 / 15</p>
---	---

Fabrication Nº:	721021-70-1 / 721021-70-3	Client:	UNIPHARMA
Equipment:	VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400		
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL		MDI VCR (v1.0)	

EN

P&ID – Cooling and Heating


	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

7 ELECTRICAL AND CONTROL SYSTEM

7.1 ELECTRIC SWITCHBOARD PANEL

All electrical and pneumatic elements are located in a stainless steel electric switchboard panel, sized 1800 x 1200 x 400 mm (height x width x depth).

Outside Electric Switchboard Panel Elements	
1	Blue indicator light – ‘Manoeuvre Voltage’ indicator
2	Green indicator light – ‘In Process’ indicator
3	Red indicator light – ‘General Alarm’ indicator
4	Acustic Alarm
5	Emergency Stop Button
6	Panel Voltage Activation switch

Below are described the exterior elements of the electric switchboard panel:

1 – BLUE INDICATOR LIGHT (MANOUEVRE VOLTAGE)

The indicator will remain lit if there is no anomaly in electrical supply.

2 – GREEN INDICATOR LIGHT ('IN PROCESS')

The indicator will remain lit if there is no anomaly in electrical supply.

3 – RED INDICATOR LIGHT ('GENERAL ALARM')

The indicator will blink when there are active alarms.

4 – ACUSTIC ALARM

Emits an audible signal whenever any alarm is activated.

5 – EMERGENCY STOP BUTTON

When pressed, stops automatically the equipment.

6- PANEL VOLTAGE ACTIVATION SWITCH

Allows the activation of the electrical switchboard panel for operation.

7.2 ELEMENTS INSIDE PANEL

Below there is the description of the most important elements located inside the electric switchboard.

FREQUENCY INVERTER

There are four frequency inverter to control the speed of VCR and MET mixers.

	<u>Anchor mixer (VCR)</u>	<u>Anchor mixer (MET)</u>
Manufacturer:	ABB	ABB
Model:	ACH550-01-08A8-4	ACH550-01-04A1-4
Tension:	380 - 480 V	380 - 480 V
Consumption:	8,8 A	4,1 A
Power:	4 kW	1,5 kW

	<u>Homogenizer (VCR)</u>
Manufacturer:	ABB
Model:	ACH550-01-045A-4
Tension:	380 - 480 V
Consumption:	44 A
Power:	22 kW

PLC


Brand:	SIEMENS
Model:	Simatic S7-1200

For extended detail of all the electric components in the switchboard, see the list attached to the electric diagrams.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

7.3 CONTROL PANEL

The equipment has a control panel for its location inside the equipment’s room. This control panel integrates the alarm indicators, the touch control panel and an emergency stop button. It consists on a stainless steel panel with dimensions of 500x500 mm.

Exterior Control Panel Elements	
1	Red indicator light – ‘General Alarm’
2	Audible Alarm
3	Touch screen panel
4	Emergency Stop Button

Below are described the exterior elements of the control panel:

1 – RED INDICATOR LIGHT (‘GENERAL ALARM’)

	The indicator will blink when there are active alarms.
---	--

2 – AUDIBLE ALARM

	Emits an audible signal whenever any alarm is activated.
---	--

3 – TOUCH SCREEN PANEL

	<p>It allows to visualize and control the status of the equipment, the alarms and also to edit the configurable parameters.</p> <p>Brand: Siemens Model: HMI TP1200 Comfort</p>
---	---

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

4 – EMERGENCY STOP BUTTON

	When pressed, stops automatically the equipment.
---	--

7.4 COMMUNICATION SYSTEM

Communication between PLC and actuators is performed by the communication system AS-interface (AS-i). This system is designed for connecting simple field I/O devices in discrete manufacturing and process applications using a single 2-conductor cable.

Otherwise, all valve actuators return a feedback of position through the same AS-i bus, in order to control its proper operation.

	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

8 INSTALLATION

8.1 SECURITY WARNINGS

Installation, start-up and maintenance may be performed by properly trained and qualified staff.

During the installation, commissioning and maintenance it is necessary to use the personal protective equipment.

It is essential to follow all safety rules when working, understand the precautions listed in technical documentation, as well as any other relevant security precaution.

Electrical connections should be done according to the electrical regulations in force and should only be handled by qualified staff.

8.2 TRANSPORT

For equipment transport, the following warnings must be taken into account:

Do not turn the packaging box upside down.

The equipment must always be transported in horizontal position.

Material should be stored in a dry, fresh place away from light.

8.3 EQUIPMENT INSTALLATION

8.3.1 Previous verification

Before installing and connecting the equipment, the following verification must be carried out:

Verify that the network's voltage matches the voltage indicated on the equipment plate.

Verify that the service connections necessary to start-up the equipment are available, and that they have appropriate characteristics.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

8.3.2 Equipment installation

To install the equipment, proceed as follows:

1. First of all, place the equipment in vertical position in its final location.
2. Install weighing modules in VCR vessel, according to instructions detailed in section 8.3.3.
3. Install the other elements supplied separately (venting filter, rupture disc, spray balls,...) in its own location. Verify that all elements are installed in its proper location according to the drawings 008601-E-0411-024-000 and 008601-E-0411-025-000.
4. Place the homogenizer on the ground.
5. Make connections between homogenizer and VCR vessel.
6. Make connections between VCR vessel, MET tank and skid.
7. Connect the equipment to the electric switchboard panel and panel control.
8. Connect the electric switchboard panel to electric supply.

8.3.3 Installation of Weigh Modules in VCR 1002 vessel

The equipment is supplied with the weigh modules disassembled in order to prevent damages during transport.

For installation of weigh modules, proceed as follows (see figure 4):

1. Raise the vessel and place the weigh modules in the vessel's legs as shown in figure 5.
2. Lower the vessel to the floor, with special care.
3. Loosen the B bolt until you can remove the centering washer
4. Tighten the B bolt.
5. Take out the A bolt.
6. Raise the vessel slowly and very carefully until the blocking plate is unlocked and you can remove it.
7. Lower the vessel carefully.

Weigh modules should be connected according to figure 5, **taking care of the cells and connection wire direction**. A wrong connection of load cells can lead to erroneous weight readings.


Fig 4. View of wheigh module


Fig 5. Connection and direction of weigh modules

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9 START-UP AND OPERATION

IMPORTANT INFORMATION!

Some elements of vessel are supplied transported without installing. Once the equipment is located in its final location and BEFORE its start-up you should ensure that all elements have been installed in their own location.

The pressure regulation valve in skid for emptying the dimple jacket should be set up to a **maximum of 0'5 bar**. A greater pressure can damage the equipment.

9.1 START OF OPERATION

To start the equipment operation, it is necessary to connect it to the electrical network through the electric switchboard panel (according to plate characteristics) and turn on the main switch located in control panel.

On starting the reactor, the “voltage drop” alarm will be triggered. You can reset this alarm via the “ALARM” screen, which you can get to from the reactor’s start screen.


9.2 RECIRCULATION OPERATIONS

The vessel has a recirculation pipe with 4 manually operated arc valves which allow the vessel’s operating mode to be selected according to their position:

- Circulation with upper shearing (Fig 6.1)
- Circulation without upper shearing (Fig 6.2)
- Circulation with lower shearing (Fig 6.3)
- Circulation without lower shearing (Fig 6.4)
- Circulation with discharge shearing (Fig 6.5)
- Circulation without discharge shearing (Fig 6.6)
- Circulation with CIP shearing (Fig 6.7)
- Circulation without CIP shearing (Fig 6.8)

In all those operations requiring shearing, make sure the manually operated valve at the back of the homogeniser is open.

Below there are figures describing the positions of the manually operated valves for each operation. The path enable with each operation is marked in blue:


Fig 6.1. Circulation with upper shearing


Fig 6.2. Circulation without upper shearing


Fig 6.3. Circulation with lower shearing


Fig 6.4. Circulation without lower shearing


Fig 6.5. Circulation with discharge shearing


Fig 6.6. Circulation without discharge shearing


Fig 6.7. Circulation with CIP shearing


Fig 6.8. Circulation without CIP shearing


Once you have selected the reactor’s operating mode and before you start it, the programme will display the image corresponding to the selected operation and request you to check and confirm that the position of the valves is correct.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.3 PRELIMINARY SETTINGS

9.3.1 Setting the drain times

Before starting-up the equipment, **you should set up the drain times** to ensure the perfect emptying of the dimple jacket before inserting steam in the heating phase. Whenever the status of the vessel changes from cooling to stop or heating, the valves involved in drainage opens during the time configured in the submenu 'SETTINGS'.

To configure drain times, proceed as follows:

1. Set up drainage duration to 0 seconds.
2. Ensure the vessel is empty, verifying that the weight indicator shows 0.
3. Start the cooling process for few minutes.
4. Stop the cooling process.
5. Manually empty the dimple jacket using the button Empty in Overview screen (Temperature Ctrl), counting the time until the weight indicator remains stable to 0.
6. The counted time will be the value which has to be set up in drain times, addend some security seconds.

9.3.2 Spill setting

In order for the product to be properly dosed, first you must set the spill parameter at the SERVICE SETTINGS screen under the SERVICE tab.


The spill (in kg) is the amount of product which will remain in the dosing tube or hose once the dosing valve is closed and which will therefore also go into the reactor.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.3.3 Setting the maximum weight oscillation

The weighing system will scan the weight for 10s in order to ensure no product is going inside the reactor and the weight remains stable.

Before starting the reactor, you must set the maximum admissible oscillation value for considering the weight to be stable.

9.4 EXTERNAL MEMORY

The touchscreen has an SD memory card inserted in the back (inside the electrical panel). This memory card will store all the recipe preparation records and the alarm history.

It is advisable to periodically dump the content of the card to avoid slowing down the system.

IMPORTANT NOTE:

The system cannot run without the SD card for storing the corresponding records. If you remove the card or if it becomes full and cannot store any more records, the programme will trigger a “Record failure” alarm.

You must restart the equipment every time you insert the SD card in the touchscreen in order for it to be recognised.

9.5 LIGHT OPERATION

The light is turned on and off by directly pressing the button on the top of the projector (on the vessel).


	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.6 TOUCH-SCREEN OPERATION

Once the equipment, the switchboard electric panel and the control panel are properly connected, the system is controlled by the touch-screen located in control panel.

9.6.1 Program structure



Program control has a main screen through which the user can select the equipment that wants to control/visualize.

The main screen also offers the following options:

OPEN/CLOSE	Log in/out the actual session.
CONFIG.	Access to HMI settings' screen.
ALARMS	Access to alarms' screen.

Access to the screens of control program is performed by pressing on the buttons located in the left part of the main screen (under the label PROCESS). Once the user has acceded, appears a bottom bar which allows the navigation through program screens:

	START	Access to the Main Menu
	OVERVIEW	Access to the Overview screen for equipment control or recipes execution.
	ALARMS	Access to Alarms screen ⁴
	RECIPES	Access to the Recipes configuration screen (creation and modification).
	CHART	Access to the Chart screen.
	SETTINGS	Access to the parameters configuration screen.


⁴ When there are active alarms, this button appears in red background.

9.6.2 OVERVIEW screen (Process / Recipe)

Overview screen allows the visualization and control of equipment and also the execution of load recipes (see section 9.6.3. for loading recipes details)

The equipment has two operating modes: Process and Recipe.

The screen has an operation selection and control panel at the top. The reactor's picture will show the current state of every parameter and the values of the set points.

To display the details on the temperature control skid (position and active valves, enable circuit...), just press the skid box and it will enlarge.

VESSEL PARAMETERS AND OPERATIONS

The reactor has the following configurable parameters and operations in both the Process and the Recipe mode:

Configurable parameters	
Time	Desired operation time
SP Mixer	Speed and direction of rotation of the anchor stirrer
Homogeniser	Start/Stop homogeniser
SP Homogeniser	Speed of the homogeniser
Temperature VCR	Desired temperature inside VCR reactor
VCR Vacuum	Vacuum pump start/stop
Dosification	Product selection: PW/GCL/GCR
SP Dosification	Desired dosing volume
Temperature Fuser	Desired temperature inside MET tank
SP Mixer Fusor	Speed of the melting tank's mixer

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

Operations	
Temperature	No selection Heat Cool Drain double jacket
Dosing	No selection Dose PW Dose PGL Dose GCR

Types of end of process/change of step	
By time	It ends the process once the user-set time has run out.
By VCR Temperature By Fuser Temperature	It ends the process once the inside of the reactor reaches the temperature set by the user.
By dosage	It ends the process once the volume set by the user has been dosed.
By time or VCR temperature	It ends the process by time or by VCR temperature, whatever happens first.

Description of the selectable operations:

1. Heating:

The reactor starts to be heated by circulating steam through the double jacket until it reaches the temperature set by the user in the operation set point or the time set by the user runs out. Before steam is circulated through the jacket, it will be automatically drained for the duration set by the user (see Section 9.3.1).

2. Cooling:

The reactor begins to be cooled by circulating cold water through the double jacket until it reaches the temperature set by the user in the operation set point or the time set by the user runs out.

3. Drain double jacket:


The jacket is drained by the drain valves being opened for the duration set by the user in the synoptic (in this case, the time entered by the user in the programme settings for the drainage prior to the heating operation will not be taken into account).

4. Dose PW/PGL/GCR:

The dosing of the selected product into the reactor starts by the dosing valve being opened, according to the value entered by the user in the operation set point field and the spill value calculated by the user (see Section 9.3.2). The process will end once the entered amount has been dosed.

PROCESS MODE

Process mode allows you to control the reactor by setting the desired parameters and choosing the operation to be carried out every time.


Once you have set the desired parameters and chosen the right operation, just press the ► button to start the reactor.

IMPORTANT NOTE:
 Before heating or cooling the reactor, make sure you have properly set the drain times as per Section 9.3.1

RECIPE MODE

Recipe mode will only allow you to run recipes previously loaded into the PLC and to display the parameters during preparation.


To start preparing a recipe, just press the ► button.

In order to record the information corresponding to the preparation of each recipe, before a recipe starts being prepared, the programme will ask you for its batch number and name.

For more details on the creation, modification and loading of recipes please refer to Section 9.6.3.

9.6.3 RECIPE screen

On the RECIPE screen you will be able to load the recipes to be prepared and to create new recipes and/or modify existing recipes.

If you cannot see the names of the configurable parameters properly, enlarge the fields by pressing and dragging the line thereof.

Recipe parameters

For each recipe, the system will allow you to configure a total of 20 steps. In each step, you can configure the following parameters:

Message	The message to be displayed on screen before the step starts
Minutes	The minutes the step will last.
Seconds	The seconds the step will last (maximum value = 59).
Homogenizer operation	Selection of the desired operation.
SP Homogenizer (rpm)	The rotational speed of the homogeniser during the step.
SP Mixing VCR	The rotational speed of the anchor stirrer during the step.
Mixing VCR	The direction of rotation of the anchor stirrer seen from the manhole.
Temperature Ctrl	Selection of the desired operation.
SP Temperature	The desired temperature value for the step (°C).
Ctrl Dosification	Selection of the desired product for dosification.
SP Dosage	The desired dosing value for the step (kg)
Vacuum Pump	Stop/start of the vacuum unit.
SP Fuser Mixer	The rotational speed of the anchor stirrer during the step.
SP Fuser Temperature	
Step Change	The criterion for moving to the next step.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

CREATE A NEW RECIPE:

IMPORTANT:

Only Supervisor, Maintenance and Administrator users can create and edit recipes.

To create a new recipe, proceed as follows:

1. Press button  to create a new recipe.
2. Enter the name of recipe.
3. Enter all configurable parameters
4. Once all values have been configured, save recipe by pressing 

MODIFY AN EXISTING RECIPE:

To modify an existing recipe, proceed as follows::

1. Select the recipe you want to modify from the drop-down list on the upper part of the screen.
2. Modify the desired parameters.
3. Once you have finished, save the recipe by pressing the button 

LOAD A RECIPE FOR ITS EXECUTION:

To load a recipe to PLC for its execution, proceed as follows:


1. Select the desired recipe from the drop-down list on the upper part of the screen.
2. Verify the parameters to ensure it is the desired recipe.
3. Press the button  to load the recipe to PLC.
4. Access the Overview screen and verify the name shown in Recipe field corresponds to the recipe you want to execute.

To start recipe execution, press the button  in Synoptic screen (Recipe control panel).

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.6.4 ALARMS screen

On this screen you will be able to see all outstanding alarms and warnings, the date and time when they occurred, and a brief description thereof.

Displaying alarms:

Outstanding alarms and warnings will be displayed at the top of the screen by means of the **WARNING** text. The Alarms button at the bottom of the screen will be highlighted in red.

Resetting an alarm:

To reset an alarm, you must press the RESET button at the bottom of the Alarms screen.

Once you have reset an alarm, the reset date and time will be displayed on screen.

Log

The system has an alarm and warning log where you can check the last alarms that have taken place.


List of set alarms/warnings:

See list of configured alarms and warnings in section 2.4.4

9.6.5 CHARTS screen

The system will generate a set of graphs during the preparation of each recipe showing the evolution of the following parameters:

- Mixers velocity
- Homogenizer velocity
- Weight
- Temperature

By pressing on the Zoom button in the desired chart, it enlarges so that the user can observe it in more detail. For the navigation on the graph and reading of values there user can use the following buttons:

	Freeze chart
	Move chart to its start
	Move graph to the right / to the left
	Zoom in / Zoom out chart
	Display / hide cursor for values reading
	Move cursor to the right / to the left

The value corresponding to the cursor's point is shown in a table at the bottom of screen.

Curva	Erlace...	Valor	Fecha/...
rpm Agitador	Contro...	0,000000	22/11/...


	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
		MDI VCR (v1.0)

9.6.6 PROGRAM SETTINGS screen

Settings screen allows the user to configure program parameters as well as accessing to the registers, by 5 submenus: SERVICE, VF ABB, CALIBRATION, PIDS and REGISTERS.

SERVICE

To configure the basic parameters for execution of equipment operations.

EMPTYING

Emptying time: The time it takes to drain the reactor’s double jacket for the draining operation and also before heating or cooling the reactor.

DOSIFICATION

Spillover The volume of water which remains in the tube section from the valve to the reactor inlet.

Max ocilation SCAN The maximum admissible weight oscillation for considering the weight to be stable.

LOAD FACTORY SETTINGS

Load settings It allows you to load the default values.

IMPORTANT NOTE:

There should be enough drain time to ensure the reactor’s double jacket is fully drained. See Section 9.3.1 for more information on how to set this time.

Fabrication Nº: 721021-70-1 / 721021-70-3

Client: UNIPHARMA


Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400

OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL

MDI VCR (v1.0)

VF ABB

To configure frequency converter parameters and its operation limits.


VF ABB

FREQUENCY INVERTER		Default value		
		VCR	MET	Homogenizer
Mixing MIX/MAX	Minimum and Maximum operating speed [rpm]	0 – 24	0 – 71	1680 -3500
Reduction Ratio	Reduction ratio of the geared motor	72,2	24,5	1
Motor MAX	Actual maximum speed of the motor [rpm]	1770	1740	3550
Tension:	Current working voltage [V]	-	-	-
Intensity:	Current working current [A]	-	-	-
Frecuencia:	Current working frequency [Hz]	-	-	-
Velocidad:	Rotational speed of the stirrer [rpm]	-	-	-
Consumo:	Consumption of the stirrer's motor [kWh]	-	-	-
Dir. Bus:	ID address of the Modbus	-	-	-

CALIBRATION

To configure calibration parameters of equipment instrumentation.

CALIBRATION		
Parameter	Definition	Default value
VCR/Fuser Temperature		
IN	Value read by the sensor or instrument (PLC input)	---
Output	Value converted by the PLC (PLC output)	---
CAL	Calibration offset value	---
Min	Minimum scaled value	0°C
Max	Maximum scaled value	140°C
Células de pesada (Weigh cells)		
Entrada (input)	Value read by the weigh cells (PLC input)	---
Salida (output)	Value converted by the PLC (PLC output)	---
CAL	Calibration offset value	---

Calibrating the weight


The weight can be calibrated as follows:

- By means of a user-set minimum and maximum range.
- By means of a known weight.

Additionally, it also possible to set the weight to 0 (RESET) and gauge (TARE).

PID CONTROLS

Allow the configuration of PID parameters for temperature control process.

TEMPERATURE PID CONTROL


Parameter	Description	Default value
Input	Actual temperature of vessel	---
Set Point	Desired temperatura of vessel	---
Control OUT	% actual opening of steam regulation valve	---
Proportional	Proportional parameter for PID regulation	<i>Note 1</i>
Integral	Integral parameter for PID regulation	<i>Note 1</i>
Derivation	Derivative parameter for PID regulation	<i>Note 1</i>
EN Manual	Selects Manual or Automatic steam regulation valve opening.	---
SP Manual	Allows to introduce % opening of steam regulation valve in case manual mode is selected	---
Status	Shows the status of PID control	---
Error	Show if there is any fault in PID control	---

Note 1: PID values are to be set when the reactor is started after being installed in its final location.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

RECORDS

Allow the access to SD memory card for the visualization of the following records created by automatically by system: alarms history⁵ and batch records. From this submenu user can only visualize those archives with PDF extension.


⁵ Alarms history archive in PDF format is only generated when touch panel starts and it may be not updated. For updated information about alarms history it is recommended to remove SD memory card and consult the alarms history archives in CSV format.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.6.7 HMI SETTINGS screen

Allow the configuration of screen parameters.

HMI SETTINGS

Users Settings	To create and modify users, passwords and automatic log-out time (<i>see section 9.6.8 for further details</i>)
HMI Calibration	To calibrate the screen
HMI Clean	Disables the screen for a few seconds to clean it
HMI Close	Closes de screen and accesses the operative system
Acoustic signal	Enable/disable the acoustic signal of touch
Language	To select the language of the program (Spanish or English).
Change date and time	To update the date and time shown on the screen

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.6.8 USER SETTINGS screen

From User Settings screen it is possible to modify the password, automatic log-out time allows and permissions group, as well as create new users.

The system has 4 default users: Operator, Operator2, Maintenance and Admin, with the permissions described below:

Menu	Function	Username			
		Operator	Operator2	Maintenance	Admin
OVERVIEW	Display	YES	YES	YES	YES
	Start recipe	YES	YES	YES	YES
	Stop recipe	NO	YES	YES	YES
	Process operations	NO	YES	YES	YES
RECIPE	Load recipe	YES	YES	YES	YES
	Edit recipe	NO	YES	YES	YES
ALARMS	Display	YES	YES	YES	YES
	Reset	YES	YES	YES	YES
CHARTS	Display	YES	YES	YES	YES
PROGRAM SETTINGS	Edit	NO	NO	YES	YES
HMI SETTINGS	User	NO	NO	NO	YES
	Clean screen	NO	YES	YES	YES
	Calibrate screen	NO	YES	YES	YES
	Close screen	NO	NO	YES	YES
	Sound signal	NO	NO	YES	YES
	Date and time	NO	NO	YES	YES

The default passwords for each user are detailed in the following table:

User level	Username	Password
Operator	Operator	5555
Advanced operator	Operator2	4444
Maintenance	Maintenance	3333
Administrator	Admin	2222

Once the session is started, the user logged-id appears in the left upper corner in all screens.

Each user has a default automatic log-out time of 5 minutes, after which, if user has not taken any action or no buttons are pressed, the session expires and it is necessary to log-in again.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.7 COMPLIANCE WITH 21 CFR, PART 11

9.7.1 Introduction

Equipment control system has been developed according to the ‘Code of Federal Regulations’ CFR (Title 21, part 11), in compliance with the following requirements:

- Users are identified and authenticated before allowing use of the equipment.
- All actions of the users involved in equipment use for manufacturing are documented and confirmed by electronic signature.
- All operator actions requiring mandatory traceability are recorded and stored in a file (Audit Trail)

Thus it is guaranteed the complete traceability of the user who has performed each operation, the time of execution and, if necessary, the reason for the action.

Also, the system ensures that the user’s operation records are not counterfeit and can be consulted when necessary.

Equipment has the appropriate modules to comply with 21 CFR (part 11). Nevertheless, it is responsibility of the end customer to have the corresponding **SOP’s (Standard Operational Procedures)** to ensure the proper use of the system (making backups, users trainings,...) and the procedure for users managing (how to proceed in case of theft or loss of password,...)

9.7.2 Identification and authentication of users

To access the program it is required that the user is identified with a user name and password. The program comes with 4 default users with different permission levels: Operator, Operator2, Maintenance and Admin (*see section 9.6.8 for further details on user’s permissions*).

IMPORTANT!

Passwords are **personal and not** transferable. It is user’s responsibility to ensure the proper use of them.

After 3 attempts to access with wrong passwords, the user will be automatically transferred to the group of “unauthorized” and will not be permitted to log-in again.

Passwords are valid for 90 days and after, the system will require its change. In order to ensure the confidentiality of passwords, they may not be reused until have passed 270 days from its use (two consecutive password changes).

9.7.3 Automatic log-off

To avoid program access and handling by unauthorized users, each user has a configurable time for automatic log-off, after which, if user has not taken any action or no buttons are pressed, the session expires and it is necessary to log-in again.

Default log-off time is 5 minutes. It can be only modified by Administrator (ADMIN) user through User Settings screen (*see section 9.6.8*).

	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.7.4 Actions register

All operations considered as critical (those that can affect the final product) are recorded in the SD memory card located on control screen, 'Audit Trail 0' file.

Depending on the critically level of the action, prior to its execution, it is required the user to enter its electronic signature (password) and a brief comment detailing the reason for the order (see listings attached below).

Actions recorded without electronic signature requirements

- Log-in/Log off user
- Change of Process MODE / Recipe MODE
- Start/Stop Process
- Parameters configuration in Overview screen (Process MODE)
- Start/Stop Recipe
- Message confirmation
- Reset alarms
- Silence alarms
- Weight calibration. Known weight
- Weight calibration. Tare Mode
- Weight calibration. Reset Mode.

Actions recorded with electronic signature requirements

- Create/Save/Load new recipes
- Users settings
- Change date and hour
- Load Factory default values
- Parameters configuration:
 - Parameters in SERVICE screen
 - Parameters in PID screen
 - Reduction parameter in VF ABB screen
 - Temperature parameters (Min, Max, CAL) in CALIBRATION screen

Audit Trail files are important documents, so it is necessary to perform a periodically backup, by saving the files on an external PC.

it is responsibility of the end customer to have the corresponding **SOP's (Standard Operational Procedures)** to ensure the proper use of the system (making backups, users trainings,...) and the procedure for users managing (how to proceed in case of theft or loss of password,...)


	Fabrication N°: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

9.7.5 Audit Trail evaluation

Audit Trail features a protection system that ensures that the file has not been tampered with. To assess the validity of the file it is required the application **Simatic WinCC Audit Viewer** installed on an external PC.

To verify the validity of Audit Trail file, proceed as follows:

1. Install “Simatic WinCC Audit Viewer” in an external PC.
2. Remove SD memory card from control screen and insert to external PC to verify content authentication.
3. Copy the file ‘Audit Trail’ in the PC.
4. Start the application “Simatic WinCC Audit Viewer”
5. Load the file ‘Audit Trail’ by selecting in menu ‘File’ - ‘Select Database’ or by clicking directly the button Open files (see images).
6. In window, select ‘**Audit Trails/Backups/ Archives/Data Records**’. Then, search the Audit Trail file in the folder where it has been saved. Once it is located, select the file and press ‘**Open**’.

7. Automatically the file will open and the indicator ‘Data Validity’ will appear in green color if the Audit Trail file has not been tampered or in red color if any kind of manipulation has been detected. In that case, it will appear also a warning in screen.

If case of detect a manipulated Audit Trail file (Data Validity in red), it will not be valid and so, it is recommended to delete it from the SD memory card. Thus the program will create a new Audit Trail file when the SD memory card is inserted again on the screen.

it is responsibility of the end customer to have the corresponding **SOP's (Standard Operational Procedures)** to ensure the proper use of the Simatic WinCC Audit Viewer application.

Fabrication Nº: 721021-70-1 / 721021-70-3

Client: UNIPHARMA


Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400

OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL


MDI VCR (v1.0)

EN

View of an Audit Trail file without manipulation
Data Validity indicator appears in green.

View of an Audit Trail file manipulated
Data Validity indicator appears in red.

Fabrication Nº:	721021-70-1 / 721021-70-3	Client:	UNIPHARMA
Equipment:	VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400		
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL			MDI VCR (v1.0)

10 MAINTENANCE AND CLEANING

10.1 BASIC MAINTENANCE

VESSEL			
Component	Part	Operation	Periodicity
Anchor mixer	Overall	In-depth inspection: leak control, warping, corrosion... Control of tightening torques and/or loosening prevention systems	Every 12 months
	Bearing	Visual inspection of surface Replacement if necessary	Every 12 months
	Bearing seal	Replacement of seal	Every 12 months
Homogeniser	Overall	General inspection	Every 12 months
	Slip joint	Verification of the good condition of the joint Replacement if necessary	Every 6 months
	Bearings	Lubrication of bearings	Every 2500 hours
	V-belt	Check of tension	Every month
	Barrier pressure installation	Inspection of pressure gauge and check of liquid level	Every day
Manhole	Manhole seal	Visual inspection	Every week
Arc valves	Overall	Visual inspection Verification of absence of leaks and smooth operation of valve	Every month
	EPDM seals	Replacement of seal	Depending on use
Diaphragm valves	Overall	Visual inspection Verification of absence of leaks and smooth operation of valve	Every month
	Diaphragm	Replacement of diaphragm	Depending on use
Butterfly valve	Overall	Visual inspection Verification of absence of leaks and smooth operation of valve	Every month
	EPDM seals	Replacement of seal	Every 12 months
Vent filter	Filter	Visual inspection Replacement when necessary	Every month
Rupture disc	Rupture disc	Visual inspection Replacement of rupture disc	Every 12 months
Safety valve	Overall	Verification of correct operation	Every 12 months

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

TEMPERATURE CONTROL SKID			
Element	Piece	Operation	Frequency
Skid	Connections	Visual inspection.	Monthly
	Insulation	Visual inspection.	Monthly
Valves	SEE SUPPLIER MANUAL		

In any case, manufacturer manuals for each element must be consulted for more detailed instructions on use, maintenance and cleanliness.

NO repair should be done on the mixing system while there is voltage in its motor.

Before disassembling any part of the equipment, electrical supply must be disconnected from the control panel.

Before any operation/action equipment surface and its components should be cooled to avoid burns caused by direct contact with hot surfaces.

In those cases where you work with food regulations (FDA or EHEDG) and whenever there is any risk of product contamination in case of leakage, food oils and greases should be used to lubricate elements.

The indicated intervals are guidelines, since they depend on several factors. It is the client's responsibility to revise the periods to carry out each operation and to establish a maintenance plan according to the equipment's use, the products used, the working temperature, etc...

If you have any questions about the equipment's operation, contact AIRPLAN S.A. before carrying out any operation.

10.2 MOTORS

Tracking

You should check regularly that the motor is in good conditions:

- Clean enough to facilitate heat dissipation and reduce the possibility of overheating
- All parts are properly assembled (fan protection, terminals cover, cable glands, etc.)
- Motor condition is good (no bumps, rust, embedded dirt, presence of water or moisture, etc.)
- Motor works properly (no noise, excessive temperature, excessive vibration, etc.)

Bearing lubrication

Usually motor bearings are greased for life so it covers the whole life of bearings, requiring no special mechanic maintenance.

In case bearings are not greased for life there will be a regreasing device and it will be indicated in motor's plate.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
		MDI VCR (v1.0)

10.3 GEAR UNITS

Tracking

It should be checked periodically the good condition of the gear unit (no oil leaks, good condition of plugs, peepholes, no corrosion in housing, etc,...).

Lubrication

- **Gear units with lifetime lubrication:**

For normal conditions of load, humidity, temperature, cleanliness, etc. gear unit doesn't need any change of oil. Used lubricants, in normal operating conditions, reach from 20.000 to 25.000 hours of operation.

However, even if the gear unit is lifetime lubricated, in those cases of gears that may be exposed to dirt, pollution, humidity... it is recommended to perform a preventive oil change, for example, every 3 years.

- **Gear units with oil lubrication:**

In gear units with oil lubrication it should be checked periodically that oil level is good, adding oil if necessary. In any case, oils with different quality should not be mixed. In case of doubt it is recommended a complete replacement of oil.

Used grease is indicated in gear characteristics plate (*see supplier's manual for further details of brand and type of grease*).

Usual change of oil

In normal operating and temperature conditions, the following indicative frequencies can be followed:

- Mineral oil → every 4.000 hours of operation or an annual replacement (the firstly reached value).
- Synthetic oil → every 15.000 hours of operation or 3 years, (the firstly reached value)⁶.

10.4 TIGHTNESS SEALS

Inspection and replacement

All tightness seals (in manhole, covers, clamps, nozzles, sight glasses,...) should be checked periodically in order to verify that there are in perfect conditions. If any problem of conservation is detected, proceed to the seal replacement with a new one of the same material.

Assembly of a tight seal

In those tight seals pressed by screws, tightening should be done in cross for a proper and uniform fit in the whole length of the seal.

Assembly of a glued seal

If the seal is bonded to avoid that it can move or knock out of position, when removing the older seal, adhesive remains in support surface should be cleaned properly, to ensure that the new seal can be perfectly assembled.

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	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

For a proper seal bonding, enough adhesive curing time should be ensured.

Adhesive strenght increases during the first days as it cures and for this process it is very important to have a good contact wiht air.

When seals are subject to considerable efforts it is desirable to achieve a good bonding, so it should pass 24 hours after the bonding in order to reach a value close to the máximum possible.

10.5 DRY MECHANICAL SEAL

IMPORTANT!

If it is detected any loss in mechanical seal, the equipment should be stopped immediately and perform a detailed inspection of mechanical seal.

In all operations on the mechanical seal it is recommended to replace the seals, as a preventative measure against a possible damage during its manipulation or against aging or wear.

Problems and solutions of the mechanical seal

Problem	Possible solution
Seats are scratched.	Proceed lapping the seats to recover them.
Tightness seals are defective or dirty.	Proceed to replace the seals for new ones.
Another part is defective.	Try to repair it or proceed to its replacement..

Assembly of the mechanical seal

If mechanical seal is removed, when it is assembled again the following verifications should be done:

- Verify that all seals have been assembled and that they and their support surfaces are perfectly cleaned to obtain a perfect tightness.
- Sealing seats should be perfectly cleaned and in perfect conditions.
- **VERY IMPORTANT:** Seal should be assembled with a proper tightening. It should be the minimum necessary to ensure the sealing and should not be excessive to avoid shortening the useful life by excessive wear of seats or breaking by overheating. So, the tightening distances should not be modified.
- After the assembly, if it is possible, verify manually that the axis of mechanical seal rotates perfectly.

For the assembly of a dry mechanical seal you can place a little of lubrication oil in the seats (never too much) in order to give a small lubrication.

You should NOT use very viscous greases or oils, as the formed layer can affect the proper seat of the mechanical seal.

	Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
	Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
	OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	

10.6 EQUIPMENT CLEANING

Depending on the products and customer needs, may be necessary to perform a first cleaning of equipment and installation before start-up, in order to remove dirt that may have accumulated during the shipping.

Customer experience is the best method for choosing the cleaning system, products to use, frequency and intensity of cleaning to perform.

Before carrying out any cleaning operation, it is obligatory to ensure that the cleaning agents are compatible with the equipment's materials and its accessories (joints...). In case of doubt, it is recommended to remove all accessories before cleaning.

Equipment can be inspected from the sight glass or manhole. If for any reason the user decides to access to the interior of equipment through the manhole, it is **VERY IMPORTANT** to consider the following:

- Disconnect the equipment and take all necessary measures in order to prevent an accidental start of equipment while people are inside.
- Customer should know the products that use and, in case of toxic or hazardous products, should proceed to take appropriate safety measures.

Initial cleaning	
Procedure	<p>Eliminate all leftovers from polishing, cleaning the reactor's entire internal surface thoroughly with a rag.</p> <p>To eliminate stains (with a rust-like appearance) that form due to the residual humidity from tests (FAT) carried out with water, apply a damp rag, a nitric-acid solution at 3% v/v maximum to them.</p> <p>Put hot potable water (60-80º) and detergent into the reactor.</p> <p>Circulate the hot water inside the reactor, and if possible under pressure, until all waste is eliminated, as well as impurities from the detergent.</p> <p>Rinse two or three times with distilled water or demineralised water.</p>
Frequency	<p>First cleaning after installation. Mandatory for all new equipment constructed with stainless steel, before carrying out the corresponding validation process.</p>

Fabrication Nº: 721021-70-1 / 721021-70-3	Client: UNIPHARMA
Equipment: VESSEL FOR CREAMS – Mod. VCR 1002 / MELT TANK – Mod. MET 0400	
OPERATING INSTRUCTIONS AND MAINTENANCE MANUAL	
MDI VCR (v1.0)	

Radical cleaning


Procedure	<p>Prewash with hot purified or potable water (60-80°C) immediately afterward to end the manufacturing operations to avoid the later appearance of deposits and incrustations.</p> <p>Wash in open circuit with purified water, a sodium hydroxide solution or alkaline detergent, at (60-80°C), according to the type of waste to be eliminated.</p> <p>Wash in closed circuit (with recirculation) under the same conditions as the previous point.</p> <p>Rinse in open circuit with hot purified water if the sodium hydroxide or alkaline detergent solution was used.</p> <p>Wash cold under cold circuit with disinfectant solution and later rinse with purified water to eliminate disinfectant leftovers.</p> <p>Final rinse with hot purified water (60-80°C). The line must be leak-tight sealed until use.</p>
Frequency	Whenever there is a change in product, and when it is the same product and it has been established as such (weekly, after maintenance operations, etc.).

Ordinary cleaning


Procedure	<p>Prewash with hot purified or potable water (60-80°C).</p> <p>Wash in closed circuit with purified water.</p> <p>Wash in closed circuit with diluted sanitising solution. For example, sodium hypochlorite solution with 5 ppm of free chloride.</p> <p>Rinse with hot purified water immediately before starting manufacturing with the following batch.</p>
Frequency	On a daily basis, given that the product does not change.

Cleaning after period of inactivity



Procedure	<p>Wash in closed circuit with diluted sanitising solution. For example, sodium hypochlorite solution with 5 ppm of free chloride.</p> <p>Final rinse with hot purified water (70-80°C).</p>
Frequency	<p>It is considered that the cleanliness status of a piece of equipment after a radical cleaning is not indefinitely maintained. Rather, it has a validity period of from 24 hours to several weeks.</p> <p>After this period, this cleaning will be applied if it is not deemed necessary to use radical cleaning.</p>

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